

User Manual

Version 1.3 – UMIT2300



Evercode™ TCR

For use with

ECIT2300 (Human)

ECIT2301 (Human)

ECIT1310 (Mouse)

ECIT1311 (Mouse)

Legal Notices

This document and its contents are proprietary to Parse Biosciences, Inc. ("Parse Biosciences") and are intended solely for use by its customers in connection with the use of the product(s) described herein and for no other purpose. The products may be used solely FOR RESEARCH PURPOSES, AND MAY NOT BE USED IN ANY DIAGNOSTIC OR THERAPEUTIC USE IN HUMANS OR ANIMALS. This document and its contents shall not be used or distributed for any other purpose and/ or otherwise communicated, disclosed or reproduced in any way whatsoever without the prior written consent of Parse Biosciences.

No rights are granted under this document with respect to any of Parse Biosciences' intellectual property rights. The license to use of any products described herein is subject to a separate written agreement between Parse Biosciences and the applicable user.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. Parse Biosciences shall have no liability for any direct, indirect, consequential or incidental damages arising out of any failure to use the product(s) in strict compliance with the terms herein.

This document may contain references to third-party sources of information, hardware or software, products, or services and/ or third-party web sites (collectively "Third Party Information"). Parse Biosciences does not control, and is not responsible for, any Third Party Information. The inclusion of Third Party Information in this document does not imply endorsement by Parse Biosciences of the Third Party Information or the third party in any way.

The product(s) described in this document are provided for one-time use by the purchaser and may not be re-used, refurbished or resold. In addition, such product(s) may not be altered, changed or modified by anyone other than Parse Biosciences and its authorized agents, and Parse Biosciences will not be liable for any such alterations, changes or modifications. The product(s) described herein may be covered by one or more of the following patents:

U.S. Pat. No. 10,900,065

U.S. Pat. No. 11,168,355

U.S. Pat. No. 11,427,856

U.S. Pat. No. 11,634,751

U.S. Pat. No. 11,639,519

U.S. Pat. No. 11,680,283

Patents pending in the U.S. and other countries

PARSE and EVERCODE are trademarks of Parse Biosciences, Inc. The information made available herein is the sole property of Parse Biosciences and its licensors.

Copyright (c) 2024 Parse Biosciences, Inc. All Rights Reserved.

Table of Contents

| | |
|--|----|
| Overview | 5 |
| Workflow | 5 |
| Protocol Timing | 10 |
| Important Guidelines | 12 |
| Part List | 17 |
| User Supplied Equipment and Reagents..... | 21 |
| Section 1: In Situ Cell Barcoding..... | 25 |
| 1.1. Set up and Sample Counting | 25 |
| 1.2. Barcoding Round 1 | 27 |
| 1.3. Barcoding Round 2 | 30 |
| 1.4. Barcoding Round 3 | 34 |
| 1.5. Lysis and Sublibrary Generation | 37 |
| Section 2: cDNA Capture and Amplification..... | 40 |
| 2.1. cDNA Capture | 40 |
| 2.2. cDNA Template Switch | 44 |
| 2.3. cDNA Amplification | 46 |
| 2.4. Post-Amplification Purification..... | 49 |
| 2.5 cDNA Quantification | 51 |
| Section 3: WT Sequencing Library Preparation..... | 52 |
| 3.1. Fragmentation and End Prep | 52 |
| 3.2. Post-Fragmentation and End Prep Size Selection..... | 54 |

| | |
|---|----|
| 3.3. Adapter Ligation | 56 |
| 3.4. Post-Ligation Purification | 57 |
| 3.5. Barcoding Round 4 | 59 |
| 3.6. Post-Barcoding Round 4 Size Selection | 61 |
| 3.7. Sequencing Library Quantification | 63 |
| Section 4: TCR Sequencing Library Preparation | 65 |
| 4.1. TCR Amplification 1 | 65 |
| 4.2. Post TCR Amplification 1 Size Selection | 67 |
| 4.3. TCR Amplification 2 | 69 |
| 4.4. Post TCR Amplification 2 Size Selection | 72 |
| 4.5. TCR Sequencing Library Quantification | 74 |
| Appendices | 76 |
| Appendix A: Sublibrary Generation Table | 76 |
| Appendix B: Sequencing Information | 77 |
| Appendix B1: Unique Dual Index (UDI) Plate - WT Sequences | 79 |
| Appendix B2: Unique Dual Index (UDI) Plate- EC Sequences | 81 |
| Appendix C: Thermocycling Programs | 84 |
| Appendix D: TCR-Specific Troubleshooting | 90 |
| Appendix E: Revision History | 91 |

Overview


Workflow


Evercode split-pool combinatorial barcoding brings a simple workflow to large-scale single cell RNA-seq experiments. The Evercode TCR kit can profile up to 100,000 cells across up to 48 different biological samples or experimental conditions. Evercode fixation kits fix and permeabilize cells so they act as individual reaction compartments. This eliminates the need for dedicated microfluidics instrumentation. Through four rounds of barcoding, the transcriptome of each fixed cell is uniquely labeled. The four rounds of barcoding can yield a vast number of possible barcode combinations, which is more than sufficient to uniquely label up to 100,000 cells while avoiding doublets.


After sequencing, the Parse Biosciences Analysis Pipeline assigns reads that share the same four barcode combinations to a single cell.


Evercode TCR expands capabilities of scRNA-seq by simultaneously capturing T-cell receptor (TCR) information. A sufficient sample size is crucial to detecting and tracking rare clones. Discover up to 100,000 T cells, identify phenotypes and their paired TCR sequences.

The tables below provide a high-level overview of the workflow.










Section 4: TCR Sequencing Library Preparation




Protocol Timing

The table below provides details of the total and hands-on time required for the whole transcriptome workflow. A visual representation of the workflow is shown below the table.

| DESCRIPTION | TIME | HANDS-ON-TIME | SAFE STOPPING POINTS |
|---|----------------------|----------------------|----------------------------------|
| Section 1: <i>In Situ</i> Cell Barcoding | | | |
| 1.1 Set up and Sample Counting | Variable (30-90 min) | Variable (30-90 min) | |
| 1.2 Barcoding Round 1 | 75 min | 45 min | |
| 1.3 Barcoding Round 2 | 60 min | 45 min | |
| 1.4 Barcoding Round 3 | 30 min | 15 min | |
| 1.5 Lysis and Sublibrary Generation | 60 min | 45 min | -80°C ≤ 6 months |
| Section 2: cDNA Capture and Amplification | | | |
| 2.1 cDNA Capture | 60 min | 30 min | |
| 2.2 cDNA Template Switch | 120 min | 30 min | 4°C ≤ 18 hrs |
| 2.3 cDNA Amplification | 90 min | 15 min | 4°C ≤ 18 hrs in the thermocycler |
| 2.4 Post-Amplification Purification | 30 min | 30 min | |
| 2.5 cDNA Quantification | 30 min | 30 min | 4°C ≤ 48 hrs or -20°C ≤ 3 months |
| Section 3: Sequencing Library Preparation | | | |
| 3.1 Fragmentation and End Prep | 60 min | 30 min | |
| 3.2 Fragmentation and End Prep Size Selection | 30 min | 30 min | 4°C ≤ 18 hrs or -20°C ≤ 2 weeks |
| 3.3 Adapter Ligation | 30 min | 15 min | |
| 3.4 Post-Ligation Purification | 30 min | 30 min | |
| 3.5 Barcoding Round 4 | 45 min | 15 min | 4°C ≤ 18 hrs in the thermocycler |
| 3.6 Post-Barcoding Round 4 Size Selection | 30 min | 30 min | |

| DESCRIPTION | TIME | HANDS-ON-TIME | SAFE STOPPING POINTS |
|---|--------|---------------|----------------------|
| 3.7 Sequencing Library Quantification | 30 min | 30 min | -20°C ≤ 3 months |
| Section 4: TCR Sequencing Library Preparation | | | |
| 4.1 TCR Amplification 1 | 90 min | 30 min | 4°C ≤ 18 hrs |
| 4.2 Post TCR Amplification 1 Size Selection | 30 min | 30 min | -20°C ≤ 3 months |
| 4.3 TCR Amplification 2 | 60 min | 30 min | 4°C ≤ 18 hrs |
| 4.4 Post TCR Amplification 2 Size Selection | 30 min | 30 min | |
| 4.5 TCR Sequencing Library Quantification | 30 min | 30 min | 20°C ≤ 3 months |



Important Guidelines

The following guidelines provide additional information to obtain optimal performance. For additional questions not discussed below, please contact us at support@parsebiosciences.com. We also have a library of additional resources and videos on our support site at <https://support.parsebiosciences.com/>.

Sample Input

- This protocol begins with cells fixed with an Evercode Cell Fixation v3 kit. It is also compatible with samples fixed with Evercode Cell Fixation v2 kit.
- When working with mouse T cells, use the Evercode Cell Fixation (Mouse TCR) kit, which includes a murine RNase Inhibitor.
- Even if samples were counted before freezing, we strongly recommend counting cells again after thawing to account for any changes during storage and freeze thaw. Typically, a 5-15% decrease after thawing should be expected. These counts are used to determine how cells are loaded in the Round 1 Plate, and their accuracy is critical to recover the desired number of cells.
- When processing many fixed samples, we recommend aliquoting samples after fixation and counting the aliquots the day before using an Evercode TCR kits. Because aliquots have experienced comparable storage durations and freeze-thaw cycle, cell counts derived from these aliquots will provide a more accurate representation compared to counts obtained immediately following fixation. Aliquots should be thawed in a water bath set to 37°C in sets of 2-4 and counted with a hemocytometer or alternative counting device. Counts should be recorded in the Sample Loading Table, and any remaining counting aliquot should be discarded.
- Once fixed samples have been thawed, they should not be refrozen.

Cell Counting and Quality Assessment

- We recommend a hemocytometer for counting, but alternative counting devices can also be used. If possible, validate counts from alternative devices to a hemocytometer when using Evercode TCR kits for the first time.
- When first using Evercode TCR kits, we suggest saving images at each counting step.
- To assess sample quality, we recommend using viability stains like trypan blue or acridine orange and propidium iodide (AO/PI).

- Examples of trypan blue stained fixed cells are shown below. High quality fixed samples have single distinct cells with <5% cell aggregation and no debris. Higher levels of aggregation will lead to elevated doublets after sequencing. When quantifying fixed samples, it is critical to avoid counting cell debris to avoid overestimating the number of cells.

High Quality Sample



Aggregation



Debris




Figure 1: Example trypan blue stained fixed cells.

Avoiding RNase Contamination

- Take standard precautions to avoid introducing RNases into samples or reagents throughout the workflow. Always wear proper laboratory gloves and use aseptic technique.
- Although RNases are not inactivated by ethanol or isopropanol, they are inactivated by products such as RNaseZap RNase Decontamination Solution (Thermo Fisher Scientific). These can be sprayed on benchtops and pipettes.
- Filtered pipette tips should be used to reduce RNase contamination from pipettes.

Centrifugation

- A range of centrifugation speeds and durations are provided in this protocol. Optimize centrifugation conditions for each sample type to balance retention and resuspension efficiencies.
- Use a swinging bucket rotor for all high-speed centrifugation steps in this protocol. A fixed-angle rotor will lead to substantial cell loss.

Optimizing Cell Recovery

- It is critical to thoroughly resuspend the cells after centrifugation throughout the protocol. Resuspend by slowly and repeatedly pipetting up and down until no clumps are visible. Due to cell adherence to tubes, carefully pipette up and down along the bottom and sides of tubes to minimize cell loss.

- We do not recommend wide bore pipette tips as they make it difficult to resuspend cell pellets adequately.
- Ensure that the 15 mL centrifuge tubes that will be used are polypropylene, as polystyrene tubes will lead to substantial sample loss.
- The first time using an Evercode TCR kit, we recommend retaining supernatants after each barcoding step. In the unlikely event of unexpectedly high sample loss, these supernatants can be analyzed to identify points for optimization.

Sample Loading Table

- The Parse Biosciences Evercode WT Sample Loading Table v2 (Excel spreadsheet) should be completed before starting the experiment.
- If the table isn't working as expected, ensure that Macros are enabled in the Sample Loading Table. Be sure to only edit the colored cells in the table to avoid disturbing any calculations.
- In the unlikely event that samples are not concentrated enough according to the Sample Loading Table (Excel spreadsheet), choose to either:
 - Decrease the "Max number barcoded cells" until the Sample Loading Table no longer gives an error, which will maintain the desired proportion between samples but barcode fewer cells.
 - Add 14 μ L of undiluted sample into each designated well of the Round 1 Plate, which will result in more overall barcoded cells but change the desired proportion between samples.

Cell Strainers

- A cell strainer with an appropriately sized mesh should be used throughout the protocol. Although 30-40 μ m is appropriate for many T cell types, the mesh size should be chosen based on the sample type.
- To maximize cell retention with cell strainers, press the pipette tip directly against the mesh. Ensure ample pressure is applied to hold contact between the tip and the strainer to force liquid through in ~1 second. An example video can be found in our support site.

Vortexing

- Unless specified in an individual step, we strongly discourage vortexing samples or enzymes throughout the protocol.

Plate Sealing

- While sealing or unsealing 96 well plates, do not splash liquid onto the PCR plate seal or between wells. Securing plates in PCR tube racks will minimize this occurrence. PCR plate seals may be difficult to remove. Carefully peel the PCR plate seal while applying downward pressure on the plate to keep it in the PCR tube rack.

Magnetic Racks and Bead Cleanups

- The Parse Biosciences Magnetic Rack (Parse Biosciences) uses powerful rare earth magnets for rapid and efficient magnetic bead purifications for 0.2 mL tubes. The rack has high and low magnet positions important for optimal yield at key steps. We do not recommend substituting alternative racks.
- To alternate between the positions, the rack can be flipped upside down so the magnet is closer to the top (high) or bottom (low) of the 0.2 mL tubes. See the figure below for clarification.




Figure 2: Parse Biosciences Magnetic Rack at (A) high and (B) low magnet positions.

- To ensure material is not lost during bead purifications, ensure supernatants are completely clear before moving to the next step. The incubation times at each step are recommendations, but visual confirmation of clearing should be used to make the final determination. See the figure below for an example of completely cleared supernatants. Discarding any beads in supernatants will result in a reduction of transcripts and genes detected per cell.



Figure 3: Clear supernatant with compact bead pellets.

Sublibrary Loading

- The Evercode TCR kit generates 8 sublibraries with distinct Illumina indexing barcodes, which can be processed and sequenced independently or concurrently after Section 1 of the protocol.
- The number of cells per sublibrary is determined when the cells are divided into sublibraries in Section 1.5 according to the guidelines in Appendix A.
- When working with a new sample type, it can be beneficial to process a single sublibrary to optimize PCR conditions before running the remaining sublibraries.
- Sublibraries can be loaded with different numbers of cells, and the maximum number of cells that can be analyzed is the sum of cells across all sublibraries.
- Asymmetric sublibrary loading can enable cost-effective sequencing quality control. One sublibrary can be loaded with a few hundred cells and sequenced very deeply. This data can be used to choose an appropriate sequence depth for the remaining sublibraries.

Indexing Primers

- UDI Plate-WT and UDI Plate-EC are 96-well plates containing 48 unique dual indexing (UDI) primers. Each well is a single-use reaction sufficient for one indexing PCR reaction for a single sublibrary. UDI Plate-WT should be used for Section 3 and UDI Plate-EC should be used for Section 4. Each UDI plate can be used for multiple TCR kits.
- The UDI Plates are sealed with a pierceable foil to minimize cross-contamination. The plate seal should be wiped with 70% ethanol and pierced with a new pipette tip immediately prior to use. Avoid splashing or mixing the liquid between individual wells. Once a well has been used, it should not be resealed or reused.
- We recommend selecting UDIs column-wise from left to right (starting with indices 1-8). UDI sequences can be found in Appendix B.







Thermocycling Programs

- We recommend pre-programming the thermocycling programs required for the entire workflow. A summary of these can be found in Appendix C.

Part List

The Evercode TCR kit includes the following boxes. Safety Data Sheets for these reagents can be provided upon request.

-20°C Reagents. Store -20°C, PN HTT100 or MTT100

| LABEL | ITEM | PN | FORMAT | QTY |
|---|-------------------------|---------------------|-----------------------------------|-----|
|  | Round 1 Plate | HTT101 or MTT101 | Green semi-skirted 96 well plate | 1 |
|  | Round 2 Plate | WT102 | Blue semi-skirted 96 well plate | 1 |
|  | Round 3 Plate | WT103 | Yellow semi-skirted 96 well plate | 1 |
|  | Resuspension Buffer | WT104 | 5 mL tube | 1 |
|  | Sample Dilution Buffer | WT105 | 2 mL tube | 1 |
|  | Round 2 Ligation Buffer | WT106 | 5 mL tube | 1 |
|  | Round 2 Ligation Enzyme | WT107 | 1.5 mL tube | 1 |
|  | Round 2 Stop Buffer | WT108 | 2 mL tube | 1 |
|  | Round 3 Stop Buffer | WT109 | 5 mL tube | 1 |
|  | Pre-Lysis Wash Buffer | WT110 | 5 mL tube | 1 |

| LABEL | ITEM | PN | FORMAT | QTY |
|---|---------------------------|-------|-------------|-----|
|  | Round 3 Ligation Enzyme | WT111 | 1.5 mL tube | 1 |
|  | Pre-Lysis Dilution Buffer | WT112 | 2 mL tube | 1 |
|  | Lysis Enzyme | WT113 | 1.5 mL tube | 1 |
|  | Bead Wash Buffer | WT114 | 5 mL tube | 1 |
|  | Wash Buffer 1 | WT115 | 5 mL tube | 1 |
|  | Wash Buffer 2 | WT116 | 5 mL tube | 1 |
|  | Capture Enhancer | WT117 | 1.5 mL tube | 1 |
|  | Binding Buffer | WT118 | 1.5 mL tube | 1 |
|  | Wash Buffer 3 | WT119 | 5 mL tube | 1 |
|  | Template Switch Buffer | WT120 | 1.5 mL tube | 1 |
|  | Template Switch Enzyme | WT121 | 1.5 mL tube | 1 |
|  | Template Switch Primer | WT122 | 1.5 mL tube | 1 |
|  | cDNA Amp Mix | WT123 | 1.5 mL tube | 1 |

| LABEL | ITEM | PN | FORMAT | QTY |
|--|------------------------------|---------------------|-------------|-----|
|  | cDNA Amp Primers | HTT124 or MTT124 | 1.5 mL tube | 1 |
|  | Fragm/End Prep Buffer | WT125 | 1.5.mL tube | 1 |
|  | Fragm/End Prep Enzymes | WT126 | 1.5 mL tube | 1 |
|  | Ligation Adapter | WT127 | 1.5 mL tube | 1 |
|  | Adapter Ligation Buffer | WT128 | 1.5 mL tube | 1 |
|  | Adapter Ligation Enzyme | WT129 | 1.5 mL tube | 1 |
|  | Library Amp Mix | WT130 | 1.5 mL tube | 1 |
|  | Receptor PCR Mix 1 | HTT131 or MTT131 | 1.5 mL tube | 1 |
|  | Receptor PCR Mix 2 | HTT132 or MTT132 | 1.5 mL tube | 1 |
|  OR  | HT Primer or MT Primer | HTT133 or MTT133 | 1.5 mL tube | 1 |

4°C Reagents. Store 4°C, PN HTT200 or MTT200

| LABEL | ITEM | PN | FORMAT | QTY |
|---|--------------------|-------|-------------|-----|
|  | Spin Additive | WT201 | 1.5 mL tube | 1 |
|  | Lysis Buffer | WT202 | 1.5 mL tube | 1 |
|  | Streptavidin Beads | WT203 | 1.5 mL tube | 1 |

User Supplied Equipment and Reagents

The following materials and equipment are required to perform the protocol but are not provided within the kit. Note that this list does not include standard laboratory equipment, such as freezers.



Note: UDI Plates are included in ECIT2301 or ECIT1311, but not in ECIT2300 or ECIT1310. If using ECIT2300 or ECIT1310, UDI Plates must be purchased separately.

Equipment

| ITEM | SUPPLIER | PN | NOTES |
|--|-------------------|---------|---|
| Parse Biosciences Magnetic Rack | Parse Biosciences | SB1004 | We do not recommend alternative 0.2 mL PCR tube racks, as they may result in lower transcript and gene detection. |
| Centrifuge with Swinging Bucket Rotors | Various Suppliers | Varies | Compatible with 15 mL centrifuge tubes and 96 well plates and capable of reaching 4°C. |
| Microcentrifuge | Various Suppliers | Varies | Compatible with 1.5 mL and 0.2 mL tubes. |
| Water bath | Various Suppliers | Varies | Or equivalent thermomixer, heat block, or bead bath capable of holding temperature at 37°C. |
| Hemocytometer | Sigma-Aldrich® | Z359629 | Or other cell counting device. We recommend validating alternatives relative to a hemocytometer. |
| PCR tube rack | Various Suppliers | Varies | Capable of holding semi-skirted 96 well PCR plates and a tight fitting lid. |
| Plate Seal Applicator | Various Suppliers | Varies | Capable of adhering plate sealing films to 96 well plates. |
| Single Channel Pipettes: P20, P200, P1000 12-channel: P20, P200 | Various Suppliers | Varies | Or 8-channel pipettes can be substituted for 12-channel pipettes. |

| ITEM | SUPPLIER | PN | NOTES |
|---------------------------------|---------------------------|-------------|--|
| T100 Thermal Cycler | Bio-Rad Laboratories® | 1861096 | Or an equivalent thermocycler compatible with semi-skirted 96 well plates, 0.2mL PCR tubes with up to 100 µL of sample volume, sample heating and cooling from 4-98°C, and a heated lid capable of 30-105°C. |
| 6-Tube Magnetic Separation Rack | New England Biolabs® | S1506S | Or an equivalent magnetic rack for 1.5 mL tubes. |
| Vortex-Genie 2® | Scientific Industries® | SI-0236 | Or alternative vortex and adapter for 96 well plates, or a thermomixer or alternative shaker that can be set to 800-1000 RPM. |
| 6-inch Platform | Scientific Industries | 146-6005-00 | |
| Microplate Foam Insert | Scientific Industries | 504-0235-00 | |
| Qubit™ Flex Fluorometer | Thermo Fisher Scientific® | Q33327 | Or an equivalent fluorometer. |
| 2100 Bioanalyzer | Agilent® | G2939BA | Choose one. |
| 4200 TapeStation System | Agilent | G2991BA | |

Consumables

| ITEM | SUPPLIER | PN | NOTES |
|----------------|-------------------|---------|---|
| UDI Plate - WT | Parse Biosciences | UDI1001 | Each 96 well plate contains 48 unique single-use reactions sufficient for multiple Evercode Whole Transcriptome kits (48 sublibraries). UDI Plate - WT can be purchased separately or bundled with Evercode Whole Transcriptome kits. |

| ITEM | SUPPLIER | PN | NOTES |
|---|---------------------|---|---|
| UDI Plate - EC | Parse Biosciences | UDI1002 | Each 96 well plate contains 48 unique single-use reactions sufficient for multiple TCR sequencing libraries (48 sublibraries). UDI Plate - EC can be purchased separately or bundled with the kit. |
| Reagent basins | Various Suppliers | Varies | Sterile, nuclease-free reagent basins. |
| SWiSH™ Mini Cell Strainer | Stellar Scientific® | TC70-SWM-20 TC70-SWM-40 TC70-SWM-70 TC70-SWM-100 | Choose one or an equivalent sterile cell strainer with an appropriate mesh size for the cell type(s) being fixed (20 µm, 40 µm, 70 µm, 100 µm). We do not recommend FlowMi Cell Strainers (SP Bel-Art). |
| pluriStrainer® Mini | pluriSelect® | 43-10020-40 43-10040-40 43-10070-40 43-10100-40 | |
| Falcon® Cell Strainer | Corning® | 431750 431751 431752 | |
| EASYstrainer™, small | Greiner Bio-One™ | 542120 542140 542170 542100 | |
| Falcon® High Clarity PP Centrifuge Tubes, 15 mL | Corning® | 352097 | |
| DNA LoBind® Tubes, 1.5 mL, Snap Cap | Eppendorf® | 022431021 | Or equivalent DNA low-binding, nuclease-free 1.5 mL tubes. |
| TempAssure® PCR 8-Tube Strips, 0.2 mL | USA Scientific® | 1402-4700 | Or equivalent nuclease-free 0.2 mL PCR tubes. |

| ITEM | SUPPLIER | PN | NOTES |
|--|--------------------------|---|---|
| AMPure® XP Reagent | Beckman Coulter® | A63880 (5 mL) A63881 (60 mL) | Choose one. We do not recommend substituting other magnetic beads. |
| SPRIselect Reagent | Beckman Coulter | B23317 (5mL) B23318 (60mL) | |
| KAPA® Pure Beads | Roche® | KK8000 (5 mL) KK8001 (30 mL) | |
| SealPlate® | Excel Scientific | 100-SEAL-PLT | Or equivalent PCR plate seals. |
| Pipette Tips TR LTS 20 µL, 200 µL, 1000 µL | Rainin® | 17014961 17014963 17014967 | Or appropriate DNA low-binding, DNase/RNase-free, and filtered pipette tips. Do not use wide bore tips. |
| RNaseZap™ RNase Decontamination Solution | Thermo Fisher Scientific | AM9780 | Or equivalent RNase decontamination solution. |
| Ethyl alcohol, Pure | Sigma-Aldrich | 459844 | Or equivalent 100% non-denatured ethanol. |
| Nuclease-Free Water | Sigma-Aldrich | W4502 | Or equivalent nuclease-free water. |
| Trypan Blue | Various Suppliers | Varies | Or alternative dyes to assess cell viability, such as AO/PI. |
| Qubit dsDNA HS (High Sensitivity) Assay Kit | Thermo Fisher Scientific | Q33230 (100 assays) Q33231 (500 assays) | Or equivalent DNA quantifier |
| High Sensitivity DNA Kit | Agilent | 5067-4626 | Choose one that corresponds to the chosen Bioanalyzer or TapeStation. |
| High Sensitivity D5000 ScreenTape and Reagents | Agilent | 5067-5592 (screen tape) 5067-5593 (sample buffer and ladder) | |

Section 1: In Situ Cell Barcoding

1.1. Set up and Sample Counting

Prior to barcoding, cells are thawed and counted. Appropriate dilutions, loading concentrations, and loading positions are determined by the Sample Loading Table.

To set up for barcoding:

1. Open the “Evercode WT Sample Loading Table v2” (Excel spreadsheet), which will guide sample dilutions and plate loading in later steps.
2. Cool a centrifuge with swinging bucket rotors to 4°C.
3. Set a water bath to 37°C.
4. Fill a bucket with ice.
5. Prepare a hemocytometer, flow cytometer, or other cell counting device.
6. Gather the following items and handle as indicated below.

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|---------------------------|----------------|-----|--|
| Round 1 Plate | -20°C Reagents | 1 | Place directly on ice. |
| Round 2 Plate | -20°C Reagents | 1 | |
| ○ Resuspension Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by inverting 3x |
| ● Sample Dilution Buffer | -20°C Reagents | 1 | |
| ○ Round 2 Ligation Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice until use in section 1.3. Mix by inverting 3x before use. |
| ● Spin Additive | 4°C Reagents | 1 | Keep at room temperature. |

7. Thaw the previously fixed cell samples in a water bath set to 37°C until all ice crystals dissolve. Thoroughly mix each sample by pipetting and store on ice.

- While minimizing time on ice, count the number of cells in the sample with a hemocytometer or alternative cell counting device. Record the cell count.



Note: When processing many fixed samples, we recommend aliquoting samples after fixation and counting the aliquots the day before using an Evercode TCR kit. See the Important Guidelines section for details.

- Record the sample names and cell count in the Sample Loading Table.



Note: If there are errors in the Sample Loading Table, see the Important Guidelines section for suggestions.

- Place the Round 1 Plate into a thermocycler and run the following program.

| THAW ROUND 1 PLATE | | |
|--------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 10 min | 70°C | 26 µL |
| Step | Time | Temperature |
| 1 | 10 min | 25°C |
| 2 | Hold | 4°C |

- Based on the values defined in the Sample Loading Table, dilute each sample with
 - Sample Dilution Buffer and store on ice.
- Proceed immediately to Section 1.2.

1.2. Barcoding Round 1

Samples are loaded into Round 1 Plate. An *in situ* reverse transcription reaction adds well-specific barcodes that also act as sample barcodes. Cells are then pooled, centrifuged, and resuspended.

To add round 1 barcodes:

1. Gently remove the Round 1 Plate from the thermocycler, place in a 0.2 mL tube rack, and centrifuge for **1 minute** at 100 x g at 4°C.
2. Remove the Round 1 Plate from the centrifuge, place in a PCR tube rack, remove the plate seal, and store on ice.



Note: PCR plate seals may be difficult to remove. Carefully peel the seal while applying downward pressure to keep the plate from moving to minimize cross-contamination.

3. With the Round 1 Plate on ice, add **14 µL** of each diluted sample to the appropriate wells of Round 1 Plate as defined in the Sample Loading Table. Mix immediately after dispensing each sample by pipetting 3x.



CRITICAL! Do not reuse any tips during this step. Different tips must be used when pipetting cells into each well of the plate. Never place a tip that has entered one of the wells into a sample or different well.



Note: When pipetting the same sample into many wells, the sample must be mixed by gentle pipetting prior to each transfer to avoid cells settling. Do not vortex the samples.

4. While secured in a PCR tube rack on a flat surface, add a new plate seal.


- Place the Round 1 Plate into a thermocycler and run the program below. Upon completion, proceed immediately to the next step.

| BARCODING ROUND 1 | | | |
|-------------------|-----------------|---------------|--------|
| Run Time | Lid Temperature | Sample Volume | |
| 40 min | 70°C | 40 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 10 min | 50°C | 1 |
| 2 | 12 s | 8°C | 3 |
| 3 | 45 s | 15°C | |
| 4 | 45 s | 20°C | |
| 5 | 30 s | 30°C | |
| 6 | 2 min | 42°C | |
| 7 | 3 min | 50°C | |
| 8 | 5 min | 50°C | 1 |
| 9 | Hold | 4°C | 1 |

- Remove the Round 1 Plate from the thermocycler, place it in a PCR tube rack, and store it on ice.
- Place the Round 2 Plate into a thermocycler and run the following program. Proceed immediately to the next step while the program is still running.

| THAW ROUND 2 PLATE | | |
|--------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 10 min | 70°C | 10 µL |
| Step | Time | Temperature |
| 1 | 10 min | 25°C |
| 2 | Hold | 4°C |

8. While secured in a PCR tube rack on a flat surface, remove the plate seal from the Round 1 Plate.
9. With the plate and tube on ice, pool all wells from the Round 1 plate into a 15 mL centrifuge tube as follows:



- i. With a multichannel P200 set to 30 μL , mix the sample in row B by pipetting 3x.
 - ii. Transfer **30 μL** from row B to row A.
 - iii. Repeat i-ii for rows C-D to mix the sample then transfer to row A.
 - iv. Transfer any residual liquid in rows B-D to row A with a multichannel P20 set to 10 μL .
 - v. Ensure the cells in row A are in suspension as described in i. Then, transfer the total volume of each well in row A into the same 15 mL tube with a single channel P200 set to 200 μL .
10. Add **9.6 μL** of ● Spin Additive to the 15 mL tube with pooled cells. Do not discard the ● Spin Additive as it will be needed in another step.
 11. Mix by gently inverting the tube just once.
 12. Centrifuge the 15 mL tube in a swinging bucket rotor for **5-10 minutes** at 200-500 x g at 4°C. Immediately move to the next step after centrifugation.



CRITICAL! Ideal centrifugation speed and duration should be determined empirically for each sample type to optimize retention and resuspension efficiencies. Overly aggressive centrifugation can damage cell integrity and decrease data quality. If the centrifugation speeds used during fixation gave satisfactory retention, they should be used throughout this protocol.

CRITICAL! Move quickly and handle the samples gently to avoid dislodging the pellet, which will have a significant impact on data quality.

13. Remove the supernatant until about ~40 μL of liquid remains above the pellet. Use a P1000 for the first 1 mL and then a P200 for the remaining volume.



Note: Depending on the number of cells and cell types, a pellet may or may not be visible.

14. Fully but gently resuspend the pellet in **1 mL** of **O** Resuspension Buffer.

15. Add an additional **1 mL** of **O** Resuspension Buffer for a total addition of 2 mL. Store on ice.

16. Proceed immediately to Section 1.3.



Note: If the low input fixation workflow was performed, execute the following step before proceeding to Section 1.3.

17. Pipette the sample through a cell strainer into a new 15 mL tube with a P1000. Before transfer, gently mix the cells by pipetting 2x.

1.3. Barcoding Round 2

The pooled cells are added to the Ligation Master Mix, which is loaded into the Round 2 Plate. An *in situ* ligation reaction adds a well-specific barcode to the 3' end of the cDNA. The ligation reaction is quenched with Round 2 Stop Buffer, and the cells are pooled and strained.

To add round 2 barcodes:

1. Gather the following items and handle as indicated below:

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|---------------------------|----------------|-----|---|
| ● Round 2 Ligation Enzyme | -20°C Reagents | 1 | Place directly on ice. Briefly centrifuge before use. |
| ● Round 2 Stop Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by vortexing. |
| Round 3 Plate | -20°C Reagents | 1 | Place directly on ice. |

- On ice, prepare the Round 2 Ligation Master Mix by adding the following to the samples in ○ Resuspension Buffer prepared in Section 1.2. Mix thoroughly by pipetting 10x with a P1000 set to 1000 μ L. Store on ice.



Note: Prior to use, mix by inverting 3x the ○ Round 2 Ligation Buffer that was gathered and stored on ice in Section 1.1.

| ROUND 2 LIGATION MASTER MIX | |
|-------------------------------|------------|
| Sample in Resuspension Buffer | 2 mL |
| ○ Round 2 Ligation Buffer | 1.95 mL |
| ● Round 2 Ligation Enzyme | 20 μ L |
| Total Volume | 3.97 mL |

- Remove the Round 2 Plate from the thermocycler, place in a PCR tube rack, and centrifuge for **1 minute** at 100 x g at 4°C.
- While secured in a PCR tube rack on a flat surface, remove the plate seal from the Round 2 Plate and store on ice.
- Transfer the Round 2 Ligation Master Mix to a basin with a P1000.
- With the Round 2 Plate on ice and the basin on the bench, transfer Round 2 Ligation Master Mix to each well in the Round 2 Plate as follows:
 - Mix the sample in the basin by pipetting 2x with a multichannel P200 set to 40 μ L.
 - Transfer **40 μ L** of the mix to row A of the Round 2 Plate and mix by pipetting 2x.
 - Repeat i-ii to mix the sample in the basin then transfer to rows B-H.



CRITICAL! Do not reuse any tips during this step. Different tips must be used when pipetting cells into each well of the plate. Never place a tip that has entered one of the wells back into the basin or a different well.



Note: If the volume is insufficient to transfer the last row with a multichannel, tilt the basin and transfer the remaining volume with a single channel pipette. If the volume is still insufficient to fill every well, a few can be left empty without impacting experimental results.

- While secured in a PCR tube rack on a flat surface, add a new plate seal to Round 2 Plate.

8. Place the Round 2 Plate into a thermocycler and run the following program. Upon completion, proceed immediately to the next step.

| BARCODING ROUND 2 | | |
|-------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 50°C | 50 µL |
| Step | Time | Temperature |
| 1 | 15 min | 16°C |
| 2 | Hold | 4°C |

9. Briefly vortex the ● Round 2 Stop Buffer and ensure there is no precipitate. Transfer the entire volume of this tube to a new basin with a P1000.
10. Remove the Round 2 Plate from the thermocycler, place in a PCR tube rack, remove the plate seal, and store on ice.
11. With the Round 2 Plate on ice and the basin on the bench, transfer **10 µL** of the ● Round 2 Stop Buffer to each well in the Round 2 Plate with a multichannel P20. After each transfer, mix by pipetting exactly 3x.



CRITICAL! Do not reuse any tips during this step. Different tips must be used when pipetting ● Round 2 Stop Buffer into each well of the plate. Never place a tip that has entered one of the wells back into the basin or a different well.

12. While secured in a PCR tube rack on a flat surface, add a new plate seal to Round 2 Plate.
13. Place the Round 2 Plate into a thermocycler and run the following program. Proceed to the next step while the program is still running.

| ROUND 2 STOP | | |
|--------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 5 min | 50°C | 60 µL |
| Step | Time | Temperature |
| 1 | 5 min | 16°C |
| 2 | Hold | 4°C |

14. Place the Round 3 Plate into a second thermocycler and run the program below. Proceed to the next step while the program is still running.

| THAW ROUND 3 PLATE | | |
|--------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 10 min | 70°C | 10 µL |
| Step | Time | Temperature |
| 1 | 10 min | 25°C |
| 2 | Hold | 4°C |



Note: If a second thermocycler is not available, the same thermocycler used in step 14 can subsequently be used in step 15. However, the Round 2 Plate should be stored on ice until the Thaw Round 3 Plate program is complete.

15. Immediately upon completion of the Round 2 Stop program, transfer the Round 2 Plate from the thermocycler to a PCR tube rack, remove the plate seal, and store on ice.

16. With the Round 2 Plate on ice and the basin on the bench, transfer all the liquid in the Round 2 Plate into a new basin as follows:

- i. With a multichannel P200 set to 50 µL, mix the sample in row A by pipetting 3x.
- ii. Transfer **50 µL** from row A to the basin.
- iii. Repeat i-ii for rows B-H to mix the sample then transfer to the basin.
- iv. Transfer any residual liquid in the Round 2 Plate to the basin with a multichannel P20 set to 10 µL.

17. Pipette the sample through a cell strainer into a new basin with a P1000. Before each transfer, gently mix the cells in the basin by pipetting 2x and tilt the basin to recover as much liquid as possible.



CRITICAL! Do not directly touch the mesh of the cell strainer with gloved hands. To ensure that all of the liquid passes through the strainer, press the tip of the pipette against the mesh to create a tight seal and press the pipette plunger down steadily. All the liquid should pass through the strainer in ~1 second.



Note: Bubbles may form while pooling or straining. They will not affect the quality of the experiment.

18. Proceed immediately to Section 1.4.

1.4. Barcoding Round 3

The ● Round 3 Ligation Enzyme is added to the pooled cells, which are then loaded into the Round 3 Plate. A second *in situ* ligation reaction adds a third well-specific barcode, the Illumina Truseq Read 2 sequence, and a biotin. The sample is then pooled and strained.

To add round 3 barcodes:

1. Gather the following items and handle as indicated below:

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|-----------------------------|----------------|-----|--|
| ○ Round 3 Stop Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by vortexing. |
| ○ Pre-Lysis Wash Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by inverting 3x. |
| ● Round 3 Ligation Enzyme | -20°C Reagents | 1 | Place directly on ice. Briefly centrifuge before use. |
| ● Pre-Lysis Dilution Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by inverting 3x. |
| ● Lysis Enzyme | -20°C Reagents | 1 | Place directly on ice. Briefly centrifuge before use. |
| ● Lysis Buffer | 4°C Reagents | 1 | Place in a 37°C water bath until use. |

2. Add **20 µL** of ● Round 3 Ligation Enzyme to the basin containing the strained sample. Mix by gently pipetting 20x with a P1000 set to 1000 µL.
3. Remove the Round 3 Plate from the thermocycler, place in a PCR tube rack, and centrifuge for **1 minute** at 100 x g at 4°C.
4. While secured in a PCR tube rack on a flat surface, remove the plate seal from the Round 3 Plate.
5. With the Round 3 Plate on ice and the basin on the bench, transfer **50 µL** from the basin to each well in the Round 3 Plate as follows:
 - i. Mix the sample in the basin by pipetting 2x with a multichannel P200 set to 50 µL.
 - ii. Transfer **50 µL** of the mix to row A of the Round 3 Plate and mix by pipetting 2x.

iii. Repeat i-ii to mix the sample in the basin then transfer to rows B-H.



CRITICAL! Do not reuse any tips during this step. Different tips must be used when pipetting cells into each well of the plate. Never place a tip that has entered one of the wells back into the basin or a different well.



Note: If the volume is insufficient to transfer the last row with a multichannel, tilt the basin and transfer the remaining volume with a single channel pipette. If the volume is still insufficient to fill every well, a few can be left empty without impacting experimental results.

6. While secured in a PCR tube rack on a flat surface, add a new plate seal to Round 3 Plate.
7. Place the Round 3 Plate into a thermocycler and run the following program.

| BARCODING ROUND 3 | | |
|-------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 50°C | 60 µL |
| Step | Time | Temperature |
| 1 | 15 min | 16°C |
| 2 | Hold | 4°C |

8. Briefly vortex the ○ Round 3 Stop Buffer and ensure there is no precipitate. Transfer the entire volume to a new basin with a P1000.
9. Remove the Round 3 Plate from the thermocycler, place in a PCR tube rack, remove the plate seal, and store on ice.
10. With the Round 3 Plate on ice and the basin on the bench, transfer **20 µL** of the ○ Round 3 Stop Buffer from the basin to each well in the Round 3 Plate with a multichannel P20. After each transfer, mix by pipetting exactly 3x.



CRITICAL! Do not reuse any tips during this step. Different tips must be used when pipetting ○ Round 3 Stop Buffer into each well of the plate. Never place a tip that has entered one of the wells back into the basin or a different well.

11. Without incubation, proceed immediately to the next step.
12. With the Round 3 Plate on ice and the basin on the bench, transfer all the liquid in the Round 3 Plate into a new basin as follows:

- i. With a multichannel P200 set to 70 μL , mix the sample in row A by pipetting 3x.
- ii. Transfer **70 μL** from row A to the basin.
- iii. Repeat i-ii for rows B-H to mix the sample then transfer to the basin.
- iv. Transfer any residual liquid in the Round 3 Plate to the basin with a multichannel P20 pipette set to 10 μL .

13. Pipette the sample through a cell strainer into a new 15 mL tube with a P1000. Before each transfer, gently mix the cells in the basin by pipetting 2x and tilt the basin to recover as much liquid as possible.



Note: Bubbles may form while pooling or straining. They will not affect the quality of the experiment.

14. Proceed immediately to Section 1.5.

1.5. Lysis and Sublibrary Generation

The cell pool is centrifuged, washed, and resuspended in ●Pre-Lysis Dilution Buffer. The cells are counted and divided into sublibraries based on the Sublibrary Generation Table in Appendix A. These sublibraries are lysed and stored at -80°C.

To generate and lyse sublibraries:

1. Add **70 µL** of ● Spin Additive to the 15 mL tube with the sample. Gently invert once to mix.
2. Centrifuge the 15 mL tube in a swinging bucket rotor for **5-10 minutes** at 200-500 x g at 4°C. Immediately move to the next step after centrifugation.



CRITICAL! Move quickly and handle the samples gently to avoid dislodging the pellet, which will impact data quality.

3. Remove the supernatant until about ~40 µL of liquid remains above the pellet. Use a P1000 for the first 6 mL and then a P200 for the remaining volume.



Note: Depending on the number of cells and sample type, the pellet may or may not be visible.

4. Fully but gently resuspend the pellet in **1 mL** of ○ Pre-Lysis Wash Buffer.
5. Add an additional **3 mL** of ○ Pre-Lysis Wash Buffer for a total addition of 4 mL.
6. Centrifuge the 15 mL tube in a swinging bucket rotor for **5-10 minutes** at 200-500 x g at 4°C. Immediately move to the next step after centrifugation.



CRITICAL! Move quickly and handle the samples gently to avoid dislodging the pellet, which will impact data quality.

7. Remove the supernatant until about ~40 µL of liquid remains above the pellet. Use a P1000 for the first 3 mL and then a P200 for remaining volume.
8. Fully but gently resuspend the pellet in **60 µL** of ● Pre-Lysis Dilution Buffer for a final total volume of **100 µL**. Store on ice.



Note: Do not discard ● Pre-Lysis Dilution Buffer as it will be used in another step.

- While minimizing time on ice, count the number of cells in the sample with a hemocytometer or alternative cell counting device. Record the cell count.



CRITICAL! We strongly recommend using a hemocytometer and carefully mixing the sample before collecting an aliquot for counting to ensure accurate sublibrary loading.

- Decide how to divide cells across the sublibraries. See the “Sublibrary Loading” in the Important Guidelines section for strategies. Use the Sublibrary Generation Table in Appendix A to determine the volume of sample and ● Pre-Lysis Dilution Buffer that should be added to each sublibrary.



CRITICAL! Do not add more than 12,500 cells to a sublibrary. Adding additional cells will result in an increased multiplet rate.

- Ensure the cells are in suspension by pipetting 5x with a P200 set to 75 μL prior to each transfer. Add the appropriate volume of sample to 8 different 0.2 mL PCR tubes.
- Add the appropriate volume of ● Pre-Lysis Dilution Buffer to the 0.2 mL tubes for a total volume of 25 μL .
- Prepare the Lysis Master Mix in a new 1.5 mL tube as follows. Mix by pipetting 3x with a P1000 set to 220 μL . Store at room temperature.

| LYSIS MASTER MIX | |
|------------------|-------------------|
| ● Lysis Buffer | 220 μL |
| ● Lysis Enzyme | 44 μL |
| Total Volume | 264 μL |



Note: Ensure that there is no precipitate before using the ● Lysis Buffer.



CRITICAL! Do not place Lysis Master Mix on ice, as a precipitate will form.

- Add **30 μL** of Lysis Master Mix to each 0.2 mL tube with diluted cells. Store at room temperature.
- Vortex the 0.2 mL tube(s) for **10 seconds**. Briefly centrifuge.

16. Place the tube(s) into a thermocycler and run the following program. If continuing to Section 2 without freezing the sample, proceed to Section 2 while the program is still running.

| CELL LYSIS | | |
|------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 80°C | 55 µL |
| Step | Time | Temperature |
| 1 | 15 min | 65°C |
| 2 | Hold | 4°C |

17. Freeze the lysate(s) at -80°C or proceed to Section 2.



Safe stopping point: Sublibrary lysates can be stored at -80°C for up to 6 months.

Section 2: cDNA Capture and Amplification

2.1. cDNA Capture

The barcoded cDNA is captured with streptavidin-coated magnetic beads and washed to remove cellular debris.

To capture the cDNA:

1. Fill an ice bucket.
2. For each lysate, prepare 400 μ L of 85% ethanol with nuclease-free water.
3. Equilibrate 80 μ L of SPRI beads per lysate to room temperature.
4. Gather the following equipment:
 - i. Magnetic rack for 1.5 mL tubes
 - ii. Parse Biosciences magnetic rack for 0.2 mL PCR tubes
 - iii. Vortex with an adapter for 96 well plates
5. Gather the following items and handle as indicated below:

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|----------------------|----------------|-----|--|
| ● Streptavidin Beads | 4°C Reagents | 1 | Keep at room temperature. |
| ○ Bead Wash Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by inverting 3x. |
| ○ Wash Buffer 1 | -20°C Reagents | 1 | |
| ○ Wash Buffer 2 | -20°C Reagents | 1 | |
| ● Capture Enhancer | -20°C Reagents | 1 | Place directly on ice. Briefly centrifuge before use. |
| ● Binding Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. Mix by inverting 3x. |

6. Remove the desired tube(s) of lysate from the thermocycler (if continuing directly from Section 1) or from storage at -80°C.

- If previously frozen, incubate the tube(s) in water bath or thermocycler at 37°C for **5 minutes**.



Note: Ensure there is no precipitation before proceeding. Otherwise, an additional 5 minutes at 37°C should be sufficient.

- Briefly centrifuge and store at room temperature.
- Briefly centrifuge ● Capture Enhancer and gently mix by pipetting 2x with a P20 set to 15 µL.
- Add **2.5 µL** of ● Capture Enhancer to each tube of lysate and mix by pipetting 5x with a P200 set to 40 µL. Briefly centrifuge.
- Incubate for **10 minutes** at room temperature. Proceed immediately to the next step during the incubation.



Note: This incubation can be extended by 5 additional minutes (up to a total of 15 minutes) without negatively impacting performance.

- Vortex ● Streptavidin Beads until fully mixed. Add the appropriate volume of ● Streptavidin Beads to a new 1.5 mL tube as follows, depending on the number of lysates being processed:

| | | |
|-------------------------|-------|--------|
| Lysates being processed | 1 | 8 |
| ● Streptavidin Beads | 44 µL | 352 µL |

- Place the tube on the magnetic rack for 1.5 mL tubes until the solution clears (~2 minutes).
- Remove and discard the supernatant.
- Remove the tube from the magnetic rack and fully resuspend the bead pellet in the appropriate volume of ○ Bead Wash Buffer as follows:

| | | |
|-------------------------|-------|--------|
| Lysates being processed | 1 | 8 |
| ○ Bead Wash Buffer | 50 µL | 400 µL |



Note: Ensure no beads are stuck to the sides of the 1.5 mL tube.

16. Place the tube on the magnetic rack for 1.5 mL tubes until the solution clears (~2 minutes).
17. Remove and discard the supernatant.
18. Repeat steps 15-17 twice for a total of 3 washes.
19. Remove the tube from the magnetic rack. Fully resuspend the pellet in the appropriate volume of ● Binding Buffer as follows and store at room temperature.

| | | |
|-------------------------|------------|-------------|
| Lysates being processed | 1 | 8 |
| ● Binding Buffer | 55 μ L | 440 μ L |

20. Add **50 μ L** of ● Streptavidin Beads in ● Binding Buffer to each tube of lysate and fully mix by pipetting 2x with a P200 set to 90 μ L.
21. Place the tube(s) into a 96 well PCR tube rack, press to secure, and ensure the caps are secured tightly. Place the lid on the rack.
22. Place the rack onto a vortex mixer with a plate adapter. Push to secure. Vortex on 20% power (~800-1000 RPM) for **30 minutes** at room temperature.



Note: To ensure the beads are being mixed sufficiently, check that the beads are in solution 10 minutes into the incubation. If settled, increase the vortex mixing speed to keep the beads in solution.

23. Remove the tube(s) from the vortex mixer.
24. Briefly vortex the tube(s) on a standard vortex adapter. Briefly centrifuge without letting beads collect at the bottom of the tube(s).
25. Place the tube(s) on the high magnet position of the Parse Biosciences Magnetic Rack, so the magnet is closer to the top of the 0.2 mL tubes. Incubate until the solution clears (~2 minutes).



CRITICAL! Ensure the supernatant is completely clear before proceeding. Discarding any beads in the supernatant will result in a reduction of transcripts and genes detected per cell. See "Magnetic Racks and Bead Cleanups" in the Important Guidelines section for representative images.

26. While still on the magnetic rack, remove and discard the supernatant.

27. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **125 μ L** \bigcirc Wash Buffer 1.
28. Incubate for **1 minute** at room temperature.
29. Return the tube(s) to the high position of the magnetic rack. Incubate until the solution clears (~2 minutes).
30. While still on the magnetic rack, remove and discard the supernatant.
31. Repeat steps 27-30 once for a total of 2 washes with \bigcirc Wash Buffer 1.
32. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **125 μ L** \bigcirc Wash Buffer 2.



Note: Save \bigcirc Wash Buffer 2 to use for optional storage before cDNA amplification.

33. Incubate for **1 minute** at room temperature.
34. Proceed immediately to Section 2.2.

2.2. cDNA Template Switch

After an additional wash, the Template Switch Master Mix is added to the captured cDNA. The template switch reaction adds a 5' adapter to the cDNA.

To perform template switch:

1. Gather the following items and handle as indicated below:

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|--------------------------|----------------|-----|---|
| ○ Wash Buffer 3 | -20°C Reagents | 1 | Thaw and store at room temperature. Mix by inverting 3x. |
| ● Template Switch Buffer | -20°C Reagents | 1 | Thaw at room temperature then place on ice. Mix by inverting 3x. Briefly centrifuge before use. |
| ● Template Switch Primer | -20°C Reagents | 1 | |
| ● Template Switch Enzyme | -20°C Reagents | 1 | Keep on ice. Briefly centrifuge before use. |



Note: Ensure that there is no precipitate in the ● Template Switch Buffer before proceeding.

2. Prepare the Template Switch Master Mix in a new 1.5 mL tube as follows, depending on the number of sublibraries being processed. Mix by pipetting 10x and store on ice.

| TEMPLATE SWITCH MASTER MIX | | |
|----------------------------|-----------|--------|
| Number of Samples | 1 | 8 |
| ● Template Switch Buffer | 101.75 µL | 814 µL |
| ● Template Switch Primer | 2.75 µL | 22 µL |
| ● Template Switch Enzyme | 5.5 µL | 44 µL |
| Total | 110 µL | 880 µL |

3. Place each tube of captured cDNA from Section 2.1 on the high position of the magnetic rack. Incubate until the solution clears (~2 minutes).
4. While still on the magnetic rack, remove and discard the supernatant.

5. While still on the magnetic rack, add **125 µL** of **○ Wash Buffer 3** to each tube.



CRITICAL! Do not discard the **○ Wash Buffer 3** as it will be used in another step.

6. Incubate for **1 minute** at room temperature.

7. While still on the magnetic rack, remove and discard the **○ Wash Buffer 3**.

8. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **100 µL** of the Template Switch Master Mix.



Note: Because the Template Switch Master Mix is viscous, it may take time to fully resuspend the beads.

9. Briefly centrifuge without letting beads collect at the bottom of the tube(s).

10. Incubate for **30 minutes** at room temperature.

11. Fully resuspend each bead pellet by mixing 5x with a P200 set to 75 µL.

12. Place the tube(s) into a thermocycler and run the following program.

| TEMPLATE SWITCH | | |
|-----------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 60 min | 70°C | 100 µL |
| Step | Time | Temperature |
| 1 | 60 min | 42°C |
| 2 | Hold | 4°C |

13. Proceed immediately to Section 2.3. Alternatively, proceed to step 14 to store samples prior to cDNA amplification.

14. Place the tube(s) on the high position of the magnetic rack. Incubate until the solution clears (~2 minutes).



Note: Beads may need to be resuspended if they have settled.

15. While still on the magnetic rack, remove and discard the supernatant.

16. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **125 μ L O Wash Buffer 2.**



Safe stopping point: Template switched cDNA can be stored at 4°C for up to 18 hours. Do not freeze.

2.3. cDNA Amplification

The captured cDNA is washed and amplified with Template Switch Primer and Illumina Truseq Read 2 specific primers.

To amplify the cDNA:

1. Gather the following items and handle as indicated below:

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|--------------------|----------------|-----|---|
| ● cDNA Amp Mix | -20°C Reagents | 1 | Thaw at room temperature then place on ice. Mix by inverting 3x. Briefly centrifuge before use. |
| ● cDNA Amp Primers | -20°C Reagents | 1 | |

2. Prepare the cDNA Amplification Master Mix in a new 1.5 mL tube as follows. Mix by pipetting 10x and store on ice.

| cDNA AMPLIFICATION MASTER MIX | | |
|-------------------------------|--------------|-------------|
| Number of Sublibraries | 1 | 8 |
| ● cDNA Amp Mix | 60.5 μ L | 484 μ L |
| ● cDNA Amp Primers | 60.5 μ L | 484 μ L |
| Total | 121 μ L | 968 μ L |

3. Place each tube of template switched cDNA from Section 2.2 on the high position of the Parse Biosciences Magnetic Rack. Incubate until the solution clears (~2 minutes).



Note: You may need to pipette mix to resuspend settled beads so they separate appropriately.

4. While still on the magnetic rack, remove and discard the supernatant.

5. While still on the magnetic rack, add **125 µL** of **OWash Buffer 3** to each tube.
6. Incubate for **1 minute** at room temperature.
7. While still on the magnetic rack, remove and discard the **OWash Buffer 3**.
8. Remove tube(s) from the magnetic rack. Fully resuspend each bead pellet with **100 µL** of the Amplification Master Mix. Store on ice.
9. Determine the number of PCR cycles required for cDNA amplification based on the table below. Although these recommendations are appropriate for many cell types, the number of cycles may need to be optimized for each sample type.

| NUMBER OF PCR CYCLES | | |
|----------------------|-----------------|-------------------|
| Cells in Sublibrary | Primary T cells | Activated T cells |
| 200-1,000 | 14 | 13 |
| 1,000-2,000 | 12 | 11 |
| 2,000-6,000 | 10 | 9 |
| 6,000-12,500 | 8 | 7 |

10. Place the tube(s) into a thermocycler and run the following program.

| cDNA AMPLIFICATION | | | |
|--------------------|-----------------|---------------|------------------------|
| Run Time | Lid Temperature | Sample Volume | |
| 50-70 min | 105°C | 100 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 sec | 98°C | 5 |
| 3* | 45 sec | 65°C* | |
| 4 | 3 min | 72°C | |
| 5 | 20 sec | 98°C | Variable, see above |
| 6* | 20 sec | 67°C* | |
| 7 | 3 min | 72°C | |
| 8 | 5 min | 72°C | 1 |
| 9 | Hold | 4°C | 1 |



CRITICAL! If processing sublibraries with different numbers of cells, they should be amplified in separate thermocyclers according to the recommendations above.



Note: Annealing steps 3* and 6* have different time and temperature settings. Ensure these are correct before starting the program



Safe stopping point: Amplified cDNA can be stored at 4°C for up to 18 hours.

2.4. Post-Amplification Purification

Amplified cDNA is purified with a 0.8x SPRI bead cleanup.

To purify the cDNA:

1. For each tube of amplified cDNA, gather 400 μL of freshly prepared 85% ethanol.
2. Gather room temperature SPRI beads (80 μL per tube of amplified cDNA).



Note: Ensure the SPRI beads have been equilibrated to room temperature for at least 30 minutes.

3. Place each tube of amplified cDNA from Section 2.3 on the high position of the Parse Biosciences Magnetic Rack. Incubate until the solution clears (~2 minutes).



Note: If beads remain in solution after 2-3 minutes, pipette 3x in the bottom of the PCR tube with a P200 set to 40 μL . Then return to the magnet and incubate until the solution clears.

4. While still on the magnetic rack, transfer **90 μL** of the supernatant containing the cDNA into a new 0.2 mL tube(s). Store at room temperature.
5. Vortex the SPRI beads until fully mixed. Add **72 μL** of SPRI beads to each tube with amplified cDNA.
6. Vortex the tube(s) for **5 seconds**. Briefly centrifuge.
7. Incubate for **5 minutes** at room temperature.
8. Place the tube(s) on the high position of the magnetic rack. Incubate until the solution clears (~2 minutes).
9. While still on the magnetic rack, remove and discard the supernatant.
10. While still on the magnetic rack, add **180 μL** of 85% ethanol to each tube.
11. Incubate for **1 minute** at room temperature.
12. While still on the magnetic rack, remove and discard the supernatant.
13. Repeat steps 10-12 once for a total of 2 washes. Remove any residual ethanol with a P20.

14. While still on the magnetic rack, air dry the SPRI beads (~2 minutes).



CRITICAL! Do not over-dry the beads, which can lead to a substantial loss in yield. "Cracking" of the beads is a sign of over-drying.

15. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **25 µL** of nuclease-free water.

16. Incubate for **10 minutes** at 37°C in a thermocycler.

17. Place the tube(s) on the low magnet position of the magnetic rack, so the magnet is closer to the bottom of the 0.2 mL tubes. Incubate until the solution clears (~2 minutes).



Note: See 'Magnetic Racks and Bead Cleanups' in the Important Guidelines section for an image of the low position.

18. While still on the magnetic rack, transfer **25 µL** of the supernatant containing the purified cDNA into new 0.2 mL tube(s). Store on ice.



Safe stopping point: Amplified cDNA can be stored at 4°C for up to 48 hours or at -20°C for up to 3 months. Otherwise, proceed immediately to Section 3 and 4.

2.5 cDNA Quantification

The concentration and size distribution of the cDNA are measured with fluorescent dyes and capillary electrophoresis.

To quantify the cDNA:

1. Measure the concentration of each tube of purified cDNA from Section 2.4 with the Qubit dsDNA HS (High Sensitivity) Assay Kit according to the manufacturer's instructions. Record the concentration(s), which will be used in Section 3.
2. Assess the size distribution of each tube of purified cDNA with a High Sensitivity DNA Kit on the Agilent Bioanalyzer System or High Sensitivity D5000 ScreenTape and Reagents on the Agilent TapeStation System according to the manufacturer's instructions.




Figure 4: Expected post-amplification sublibrary cDNA size distribution. Example trace of Human cDNA run on a TapeStation (Top) and Mouse cDNA (Bottom).

Note: Samples may need to be diluted to be within the manufacturer's recommended concentration range. Typically, between a 1:3 to 1:10 dilution is appropriate.



Note: The traces above are representative of typical TapeStation cDNA traces. The shape and prominence of the trace is dependent on cell type, sublibrary size, and amount of DNA loaded into the TapeStation. Sublibraries with minor deviations can still produce high quality data.

Section 3: WT Sequencing Library Preparation

3.1. Fragmentation and End Prep

Barcoded and amplified cDNA is fragmented, end repaired and A-tailed in a single reaction.

To prepare for fragmentation and end prep:

1. For each sublibrary, prepare 1.2 mL of 85% ethanol with nuclease-free water.
2. Equilibrate 180 μ L of SPRI beads per sublibrary to room temperature.
3. Fill an ice bucket.
4. Take out the Parse Biosciences magnetic rack for 0.2 mL PCR tubes.
5. Obtain recorded cDNA concentrations from Section 2.5.
6. Gather the following items and handle as indicated below.

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|--------------------------|----------------|-----|---|
| ● Fragm/End Prep Buffer | -20°C Reagents | 1 | Thaw at room temperature then store on ice. |
| ● Fragm/End Prep Enzymes | -20°C Reagents | 1 | Place directly on ice. Briefly centrifuge before use. |

7. Vortex the tube(s) of cDNA for **5 seconds**. Briefly centrifuge.
8. Prepare Diluted cDNA in new 0.2 mL tube(s) as follows to a final volume of 35 μ L and store on ice. Store any remaining purified cDNA at -20°C.

| DILUTED cDNA | |
|---------------------|------------|
| Purified cDNA | 10 μ L |
| Nuclease-free water | 25 μ L |
| Total Volume | 35 μ L |

9. Start the following program to ensure the thermocycler is cool prior to use.

| FRAGMENTATION AND END PREP | | |
|----------------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 40 min | 70°C | 50 µL |
| Step | Time | Temperature |
| 1 | Hold* | 4°C |
| 2 | 10 min | 32°C |
| 3 | 30 min | 65°C |
| 4 | Hold | 4°C |



Note: * This hold ensures that the thermocycler is cooled and ready for step 13.

10. Vortex the ● Fragm/End Prep Buffer for **5 seconds**. Briefly centrifuge.



Note: Confirm the ● Fragm/End Prep Buffer is fully thawed without precipitation.

11. Prepare the Fragmentation and End Prep Master Mix in a new 1.5 mL tube as follows. Mix by pipetting 10x and store on ice.

| FRAGMENTATION AND END PREP MASTER MIX | | |
|---------------------------------------|---------|--------|
| Number of Sublibraries | 1 | 8 |
| ● Fragm/End Prep Buffer | 5.5 µL | 44 µL |
| ● Fragm/End Prep Enzymes | 11 µL | 88 µL |
| Total | 16.5 µL | 132 µL |

12. Add **15 µL** of Fragmentation and End Prep Master Mix to each tube of Diluted cDNA. Mix by pipetting 10x with a P200 multichannel pipette set to 40 µL. Briefly centrifuge.

13. Place the tube(s) into the cooled thermocycler. Skip the 4°C hold in Step 1 so the thermocycler proceeds to Step 2 of the Fragmentation and End Prep program.



Note: Ensure the thermocycler has cooled to 4°C prior to adding the tubes.

14. As soon as the program reaches step 4 of the thermocycling program (4°C), store the tube(s) on ice and proceed immediately to Section 3.2.

3.2. Post-Fragmentation and End Prep Size Selection

The fragmented and end prepped DNA are size selected with a double sided SPRI cleanup.

To size select the fragmented and end prepped DNA:

1. Gather freshly prepared 85% ethanol.
2. Gather room temperature SPRI beads (~50 µL per sublibrary).



Note: Ensure the SPRI beads have been equilibrated to room temperature for at least 30 minutes.

3. Vortex the SPRI beads until fully mixed. Add **30 µL** of SPRI beads to each tube of fragmented and end prepped DNA.
4. Vortex the tube(s) for **5 seconds**. Briefly centrifuge.
5. Incubate for **5 minutes** at room temperature.
6. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).
7. While still on the magnetic rack, transfer **75 µL** of the supernatant containing the fragmented and end prepped DNA into new 0.2 mL tube(s). Discard the tube(s) with bead pellet(s).
8. Add **10 µL** of SPRI beads to each tube.
9. Vortex the tube(s) for **5 seconds**. Briefly centrifuge.
10. Incubate for **5 minutes** at room temperature.

11. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~3 minutes).



CRITICAL! Ensure the solution is completely clear before proceeding.

12. While still on the magnetic rack, remove and discard the supernatant.

13. While still on the magnetic rack, add **180 µL** of 85% ethanol to each tube.

14. Incubate for **1 minute** at room temperature.

15. While still on the magnetic rack, remove and discard the supernatant.

16. Repeat steps 13-15 once for a total of 2 washes. Remove any residual ethanol with a P20.

17. While still on the magnetic rack, air dry the SPRI beads (~30 seconds).



CRITICAL! Do not over-dry the beads, which can lead to substantial loss in yield. "Cracking" of the beads is a sign of over-drying.

18. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **50 µL** of nuclease-free water.

19. Incubate for **5 minutes** at room temperature.

20. Place the tube(s) on the high position of the magnetic rack. Incubate until the solution clears (~2 minutes).

21. While still on the magnetic rack, transfer **50 µL** of the supernatant into new 0.2 mL tube(s).



Safe stopping point: The size-selected fragmented and end prepped DNA can be stored at 4°C for up to 18 hours or at -20°C for up to 2 weeks.

3.3. Adapter Ligation

Adapters with an Illumina Truseq Read 1 sequence are ligated to the 5' end of the fragmented and end prepped DNA.

To ligate adapters:

1. Gather the following items and handle as indicated below.

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|---------------------------|----------------|------------------------------|---|
| ● Ligation Adapter | -20°C Reagents | 1 | Thaw at room temperature then place on ice. Mix by inverting 3x. Briefly centrifuge before use. |
| ● Adapter Ligation Buffer | -20°C Reagents | 1 | |
| ● Library Amp Mix | -20°C Reagents | 1 | |
| UDI Plate - WT | -20°C Reagents | 1 sealed well per sublibrary | Thaw at room temperature then place on ice. |
| ● Adapter Ligation Enzyme | -20°C Reagents | 1 | Place directly on ice. Briefly centrifuge before use. |

2. Prepare the Adapter Ligation Master Mix in a new 1.5 mL tube as follows. Mix by pipetting 10x and store on ice.

| ADAPTER LIGATION MASTER MIX | | |
|-----------------------------|----------|--------|
| Number of Sublibraries | 1 | 8 |
| Nuclease-free water | 19.25 µL | 154 µL |
| ● Adapter Ligation Buffer | 22 µL | 176 µL |
| ● Adapter Ligation Enzyme | 11 µL | 88 µL |
| ● Ligation Adapter | 2.75 µL | 22 µL |
| Total | 55 µL | 440 µL |

3. Add **50 µL** of Adapter Ligation Master Mix to each tube of purified fragmented and end prepped DNA from Section 3.2. Mix by pipetting 10x with a P200 multichannel pipette set to 80 µL. Briefly centrifuge.

- Place the tube(s) into a thermocycler and run the program below.

| ADAPTER LIGATION | | |
|------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 30°C* | 100 µL |
| Step | Time | Temperature |
| 1 | 15 min | 20°C |
| 2 | Hold | 4°C |



Note: * If the thermocycler's lid can't reach the recommended temperature of 30°C, turn the lid heating off.

- As soon as the program reaches 4°C, store the adapter ligated DNA on ice and proceed immediately to Section 3.4.

3.4. Post-Ligation Purification

Adapter ligated DNA is purified with a 0.8x SPRI bead cleanup.

To purify the ligated DNA:

- Gather freshly prepared 85% ethanol.
- Gather room temperature SPRI beads (~90 µL per sublibrary)



Note: Ensure the SPRI beads have been equilibrated to room temperature for at least 30 minutes.

- Vortex the SPRI beads until fully mixed. Add **80 µL** of SPRI beads to each tube of adapter ligated DNA from Section 3.3.
- Vortex the tube(s) for **5 seconds**. Briefly centrifuge.
- Incubate for **5 minutes** at room temperature.
- Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).


7. While still on the magnetic rack, remove and discard the supernatant.
8. While still on the magnetic rack, add **180 μ L** of 85% ethanol to each tube.
9. Incubate for **1 minute** at room temperature.
10. While still on the magnetic rack, remove and discard the supernatant.
11. Repeat steps 8-10 once for a total of 2 washes. Remove any residual ethanol with a P20.
12. While still on the magnetic rack, air dry the SPRI beads (~2 minutes).
13. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **23 μ L** of nuclease-free water.
14. Incubate for **5 minutes** at room temperature.
15. Place the tube(s) on the low magnet position of the magnetic rack, so the magnet is closer to the bottom of the 0.2 mL tubes. Incubate until the solution clears (~2 minutes).
16. While still on the magnetic rack, transfer exactly **21 μ L** of the supernatant containing the purified adapter ligated DNA into new 0.2 mL tube(s). Store on ice.
17. Proceed immediately to Section 3.5.

3.5. Barcoding Round 4

Purified adapter ligated DNA is PCR amplified with Illumina Truseq Read 1 and Read 2 primers. This indexing PCR generates sequencing libraries and adds i5/i7 UDIs that act as a fourth cell barcode.

To add round 4 barcodes:

1. Centrifuge the UDI Plate - WT at 100 x g for **1 minute**.
2. Wipe the surface of the plate with 70% ethanol and allow it to dry.
3. Orient the UDI Plate - WT with the notch on the bottom left as pictured below. For each sublibrary being processed, choose one unused well of the UDI Plate - WT and record the well position and number for each sublibrary.



4. With a multichannel P20, pierce the seal of the chosen wells of the UDI Plate - WT.
5. With a multichannel P20 and new tips, mix by pipetting 5x then immediately transfer **4 μ L** from a chosen unused well of the UDI Plate - WT to its corresponding tube of purified adapter ligated DNA from Section 3.4.



CRITICAL! Only transfer primers from 1 well of the UDI Plate - WT to 1 tube of adapter ligated DNA.

6. If any unused wells remain in the UDI Plate - WT, store the plate at -20°C . Do not reuse wells.
7. Add **25 μ L** of ● Library Amp Mix to each tube. Mix by pipetting 10x with a P200 multichannel pipette set to 25 μ L. Briefly centrifuge.

8. Determine the number of PCR cycles required for the Indexing PCR based on the amount of cDNA added to the fragmentation and end prep reaction as recorded in Section 2.5.

| NUMBER OF PCR CYCLES | |
|----------------------|------------|
| cDNA Input (ng) | PCR Cycles |
| 10-24 | 13 |
| 25-49 | 12 |
| 50-99 | 11 |
| 100-299 | 10 |
| 300-999 | 8 |
| 1,000 or more | 7 |

9. Place the tube(s) into a thermocycler and run the following program.

| INDEXING PCR | | | |
|--------------|-----------------|---------------|-------------------------|
| Run Time | Lid Temperature | Sample Volume | |
| ~30 min | 105°C | 50 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | Varies, see table above |
| 3 | 20 s | 67°C | |
| 4 | 1 min | 72°C | |
| 5 | 5 min | 72°C | 1 |
| 6 | Hold | 4°C | 1 |



CRITICAL! If processing sublibraries with different cDNA concentrations, they should be amplified in separate thermocyclers according to the recommendations above.



Safe stopping point: Prior to size selection, sequencing libraries can be stored at 4°C for up to 18 hours.

3.6. Post-Barcoding Round 4 Size Selection

The sequencing libraries are size selected with a double sided SPRI cleanup.

To size select the sequencing libraries:

1. Gather freshly prepared 85% ethanol.
2. Gather room temperature SPRI beads (~50 μL per sublibrary).



Note: Ensure the SPRI beads have been equilibrated to room temperature for at least 30 minutes.

3. Vortex the SPRI beads until fully mixed. Add **30 μL** of SPRI beads to each sequencing library tube.
4. Vortex the tube(s) for **5 seconds**. Briefly centrifuge.
5. Incubate for **5 minutes** at room temperature.
6. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).
7. While still on the magnetic rack, transfer **75 μL** of the supernatant containing the DNA into new 0.2 mL tube(s). Discard the tube(s) with bead pellet(s).
8. Add **10 μL** of SPRI beads to each tube.
9. Vortex the tube(s) for **5 seconds**. Briefly centrifuge.
10. Incubate for **5 minutes** at room temperature.
11. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~3 minutes).



CRITICAL! This may take longer due to the low volume of beads. Ensure the solution is completely clear before proceeding.

12. While still on the magnetic rack, remove and discard the supernatant.
13. While still on the magnetic rack, add **180 μL** of 85% ethanol to each tube.
14. Incubate for **1 minute** at room temperature.

15. While still on the magnetic rack, remove and discard the supernatant.
16. Repeat steps 13-15 once for a total of 2 washes. Remove any residual ethanol with a P20.
17. While still on the magnetic rack, air dry the SPRI beads (~30 seconds).
18. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **20 μ L** of nuclease-free water.
19. Incubate for **5 minutes** at room temperature.
20. Place the tube(s) on the low position of the magnetic rack. Incubate until the solution clears (~2 minutes).
21. While still on the magnetic rack, transfer the supernatant into new 0.2 mL tube(s). Store on ice.



Safe stopping point: Sequencing libraries can be stored at -20°C for up to 3 months.

3.7. Sequencing Library Quantification

The concentration and size distribution of the sequencing libraries are measured with fluorescent dyes and capillary electrophoresis.

To quantify the sequencing libraries:

1. Measure the concentration of each purified sequencing library from Section 3.6 with the Qubit dsDNA HS (High Sensitivity) Assay Kit according to the manufacturer’s instructions.
2. Assess the size distribution of each purified sequencing library with a High Sensitivity DNA Kit on the Agilent Bioanalyzer System or High Sensitivity D1000 ScreenTape and Reagents on the Agilent TapeStation System according to the manufacturer’s instructions.




Figure 5: Expected Size Distribution before Illumina Sequencing. Example trace of Human DNA (Top) and Mouse DNA (Bottom) from indexed sublibraries run on a TapeStation.

Note: Samples may need to be diluted to be within the manufacturer's recommended concentration range. Typically, between a 1:3 to 1:10 dilution is appropriate.



Note: The traces above are representative of typical TapeStation of DNA from indexed sublibraries. There should be a peak between 400-500 bp. The prominence of the trace is dependent on amount of DNA loaded into the TapeStation. Sublibraries with minor deviations can still produce high quality data.

Note: If using a Bioanalyzer, there may be an additional peak present. This typically occurs if products are overamplified, but it should not impact sequencing or data quality (assuming there is still a peak present at 400-500 bp). Do not use this additional peak when estimating amplicon size.

Section 4: TCR Sequencing Library Preparation

4.1. TCR Amplification 1

A PCR reaction amplifies a subset of the cDNA containing V(D)J segments in the CDR3 repertoire of the T cells. For the first TCR Amplification Reaction, an HT Primer is used in the Human Evercode TCR, and a MT Primer is used in the Mouse Evercode.

To amplify the V(D)J segments:

1. Prepare **1.2 mL** 85% ethanol per sublibrary with nuclease-free water.
2. Equilibrate approximately **180 µL** of SPRI beads per sublibrary to room temperature.
3. Fill an ice bucket.
4. Take out the Parse Biosciences magnetic rack for 0.2 mL PCR tubes.
5. Obtain amplified cDNA from Section 2.4.
6. Gather the following items and handle as indicated below.

| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|-----------------------|--------------------|-----|---|
| ● Receptor PCR Mix 1 | TCR -20°C Reagents | 1 | Thaw at room temperature then place on ice. Mix by inverting 3x. Briefly centrifuge before use. |
| ● HT Primer/MT Primer | TCR -20°C Reagents | 1 | |

7. Vortex the tube(s) of amplified cDNA for **2-3 seconds**. Briefly centrifuge.
8. Prepare diluted cDNA in a new 0.2 mL tube(s) as follows and store on ice. Bring the total volume up to **36 µL** using nuclease-free water. Store any remaining cDNA at -20°C.

| DILUTED cDNA | |
|---------------------|----------|
| Sublibrary cDNA | 10 ng |
| Nuclease-free water | Variable |
| Total volume | 36 µL |



Note: Accuracy of input is important at this step. If the concentration of cDNA is higher than 5 ng/µL, a dilution is recommended to ensure sufficient volume is used.

Note: If you have less than 10 ng of cDNA for a given sublibrary, add the entire amount of cDNA at this step. This will not affect the quality of your libraries. Record the amount added to each tube as subsequent PCR cycles will have to be adjusted based on cDNA concentration (See Appendix D for more details).

- Prepare the TCR Amplification Reaction Solution in a new 1.5 mL tube as follows. Mix by pipetting 10x and store on ice.

| TCR AMPLIFICATION 1 MASTERMIX | | |
|-------------------------------|--------------|---------------|
| Number of Sublibraries | 1 | 8 |
| ● Receptor PCR Mix 1 | 55 μ L | 440 μ L |
| ● HT Primer/MT Primer* | 15.4 μ L | 123.2 μ L |
| Total | 70.4 μ L | 563.2 μ L |



Note: *The initial TCR Amplification Reaction uses an HT Primer for human TCR and an MT Primer for mouse TCR.

- Add **64 μ L** of the TCR Amplification Reaction Solution to each sublibrary, bringing the total volume to **100 μ L**. Mix sublibraries 10x with a P200 pipette set to 90 μ L. Briefly centrifuge.
- Place the tube(s) into a thermocycler and run the following program.

| TCR AMPLIFICATION 1 | | | |
|---------------------|-----------------|---------------|--------|
| Run Time | Lid Temperature | Sample Volume | |
| ~50 minutes | 105°C | 100 μ L | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | 20 |
| 3 | 20 s | 65°C | |
| 4 | 20 s | 72°C | |
| 5 | 5 min | 72°C | 1 |
| 6 | Hold | 4°C | 1 |



Safe stopping point: TCR enriched cDNA can be stored at 4°C for up to 18 hours.

4.2. Post TCR Amplification 1 Size Selection

TCR enriched cDNA sublibraries are size selected with a double sided 0.5X and 0.8X SPRI cleanup.

To size select the sublibraries:

1. Gather freshly prepared 85% ethanol.
2. Gather room temperature SPRI beads (~90 μL per sublibrary).



Note: Ensure the SPRI beads have been equilibrated to room temperature for at least 30 minutes.

3. Vortex the SPRI beads until fully mixed. Add **50 μL** of SPRI beads to each tube of TCR-amplified cDNA.
4. Vortex the tube(s) for **2-3 seconds**. Briefly centrifuge.
5. Incubate for **5 minutes** at room temperature.
6. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).
7. While still on the magnetic rack, transfer **145 μL** of the supernatant containing the TCR-amplified cDNA into new 0.2 mL tube(s). Discard the tube(s) with bead pellet(s).
8. Add **30 μL** of SPRI beads to each tube.
9. Vortex the tube(s) for **2-3 seconds**. Briefly centrifuge.
10. Incubate for **5 minutes** at room temperature.
11. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).



CRITICAL! This may take longer due to the low volume of beads. Ensure all the solutions are completely clear before proceeding.

12. While still on the magnetic rack, remove and discard the supernatant.
13. While still on the magnetic rack, add **180 μL** of 85% ethanol to each tube.
14. Incubate for **1 minute** at room temperature.

15. While still on the magnetic rack, remove and discard the supernatant.
16. Repeat steps 13-15 once for a total of 2 washes. Remove any residual ethanol with a P20.
17. While still on the magnetic rack, air dry the SPRI beads (~30 seconds).



CRITICAL! Do not over-dry the beads, which can lead to substantial loss in yield. "Cracking" of the beads is a sign of over-drying.

18. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **16 μ L** of nuclease-free water.
19. Incubate for 5 minutes at room temperature.
20. Place the tube(s) on the low position of the magnetic rack. Incubate until the solution clears (~2 minutes).
21. While still on the magnetic rack, transfer **15 μ L** of the supernatant into new 0.2 mL tube(s) and store on ice. Discard the tube(s) with bead pellet(s).
22. Measure the concentration of each tube of size-selected TCR amplified cDNA with the Qubit dsDNA HS (High Sensitivity) Assay Kit according to the manufacturer's instructions.



Safe stopping point: The size-selected TCR amplified cDNA can be stored at -20°C for up to 3 months.

4.3. TCR Amplification 2

A final PCR further amplifies the TCR cDNA and adds i5/i7 UDIs that act as a fourth cell barcode.

To prepare the TCR sequencing libraries:

1. Obtain amplified and purified TCR cDNA from Section 4.2.
2. Vortex the tube(s) for **2-3 seconds**. Briefly centrifuge.
3. Prepare diluted TCR cDNA in new 0.2 mL tube(s) as follows and store on ice. Bring the total volume up to 46 μ L using nuclease-free water. Store any remaining product from Section 4.2 at -20°C .

| DILUTED TCR cDNA | |
|---------------------|------------|
| TCR cDNA | 10 ng |
| Nuclease-free water | Variable |
| Total volume | 46 μ L |

Note: Accuracy of input is important at this step. If the concentration of TCR cDNA is higher than 5 ng/ μ L, a dilution is recommended to ensure sufficient volume is used while pipetting (at least 2 μ L used).



Note: If you have less than 10 ng of TCR cDNA for a given sublibrary, add the entire amount of TCR cDNA at this step. This will not affect the quality of your libraries. Record the amount added to each tube as subsequent PCR cycles will have to be adjusted based on TCR cDNA concentration (See Appendix D for more details).

4. Gather the following items and handle as indicated below.


| ITEM | SOURCE | QTY | HANDLING AND STORAGE |
|----------------------|--|------------------------------|---|
| ● Receptor PCR Mix 2 | TCR -20°C Reagents | 1 | Thaw at room temperature then place on ice. Mix by inverting 3x. Briefly centrifuge before use. |
| UDI Plate - EC | UDI Plate - EC (-20°C) | 1 sealed well per sublibrary | Thaw at room temperature then place on ice. |

5. Centrifuge the UDI Plate - EC at 100 x g for **1 minute**.



CRITICAL! Ensure not to confuse UDI plates. Double-check the label on the plate as specific plates are used in different sections.

6. Wipe the surface of the plate with 70% ethanol and allow it to dry.
7. Orient the UDI Plate - EC with the notch on the bottom left as pictured below. For each sublibrary being processed, choose one unused well of the UDI Plate - EC and record the well position and number for each sublibrary.



8. With a multichannel P20, pierce the seal of the chosen wells of the UDI Plate - EC.
9. With a multichannel P20 and new tips, mix by pipetting 5x then immediately transfer **4 μ L** from a chosen unused well of the UDI Plate - EC to its corresponding tube of diluted TCR DNA.



CRITICAL! Only transfer primers from 1 well of the UDI Plate - EC to 1 tube of diluted TCR cDNA.

10. If any unused wells remain in the UDI Plate - EC, store the plate at -20°C . Do not reuse wells.
11. Add **50 μ L** of ● Receptor PCR Mix 2 to each tube. Mix by pipetting 10x with a P200 multichannel pipette set to 90 μ L.

12. Place the tube(s) into a thermocycler and run the following program.

| TCR AMPLIFICATION 2 | | | |
|---------------------|-----------------|---------------|--------|
| Run Time | Lid Temperature | Sample Volume | |
| ~30 minutes | 105°C | 100 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | 9 |
| 3 | 20 s | 67°C | |
| 4 | 1 min | 72°C | |
| 5 | 5 min | 72°C | 1 |
| 6 | Hold | 4°C | 1 |



Safe stopping point: TCR sequencing libraries can be stored at 4°C for up to 18 hours.

4.4. Post TCR Amplification 2 Size Selection

TCR sequencing libraries are size selected with a double sided 0.5X and 0.65X SPRI cleanup.

To size select the sublibraries:

1. Gather freshly prepared 85% ethanol.
2. Gather room temperature SPRI beads (~80 μ L per sublibrary).



Note: Ensure the SPRI beads have been equilibrated to room temperature for at least 30 minutes.

3. Vortex the SPRI beads until fully mixed. Add **50 μ L** of SPRI beads to each tube of TCR sequencing library.
4. Vortex the tube(s) for **2-3 seconds**. Briefly centrifuge.
5. Incubate for **5 minutes** at room temperature.
6. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).
7. While still on the magnetic rack, transfer **145 μ L** of the supernatant containing the TCR sequencing library into new 0.2 mL tube(s). Discard the tube(s) with bead pellet(s).
8. Add **15 μ L** of SPRI beads to each tube.
9. Vortex the tube(s) for **2-3 seconds**. Briefly centrifuge.
10. Incubate for **5 minutes** at room temperature.
11. Place the tube(s) on the high position of the magnetic rack for 0.2 mL tubes. Incubate until the solution clears (~2 minutes).



CRITICAL! This may take longer due to the low volume of beads. Ensure all the solutions are completely clear before proceeding.

12. While still on the magnetic rack, remove and discard the supernatant.
13. While still on the magnetic rack, add **180 μ L** of 85% ethanol to each tube.
14. Incubate for **1 minute** at room temperature.

15. While still on the magnetic rack, remove and discard the supernatant.
16. Repeat steps 13-15 once for a total of 2 washes. Remove any residual ethanol with a P20.
17. While still on the magnetic rack, air dry the SPRI beads (~30 seconds).



CRITICAL! Do not over-dry the beads, which can lead to substantial loss in yield. "Cracking" of the beads is a sign of over-drying.

18. Remove the tube(s) from the magnetic rack. Fully resuspend each bead pellet with **21 μ L** of nuclease-free water.
19. Incubate for **5 minutes** at room temperature.
20. Place the tube(s) on the low position of the magnetic rack. Incubate until the solution clears (~2 minutes).
21. While still on the magnetic rack, transfer **20 μ L** of the supernatant into new 0.2 mL tube(s) and store on ice. Discard the tube(s) with bead pellet(s). Product is now ready for quantification.



Safe stopping point: Sequencing libraries can be stored at -20°C for up to 3 months.

4.5. TCR Sequencing Library Quantification

The concentration and size distribution of the TCR sequencing libraries are measured with fluorescent dyes and capillary electrophoresis.

To quantify the TCR sequencing libraries:

1. Measure the concentration of each purified TCR sequencing library from Section 4.4 with the Qubit dsDNA HS (High Sensitivity) Assay Kit according to the manufacturer's instructions.
2. Assess the size distribution of each purified TCR sequencing library with a High Sensitivity DNA Kit on the Agilent Bioanalyzer System or High Sensitivity D5000 ScreenTape and Reagents on the Agilent TapeStation System according to the manufacturer's instructions. Samples may need to be diluted to be within the manufacturer's recommended concentration range.




Figure 6: Expected TapeStation size distribution of a Human (Top) and Mouse (Bottom) TCR Sublibrary made with both Primary and Activated T cells.

Note: Samples may need to be diluted to be within the manufacturer's recommended concentration range. Typically, between a 1:3 to 1:10 dilution is appropriate.



Note: The trace above is representative of a typical TapeStation trace of DNA from an indexed TCR sublibrary. Amplicon sizes should be within 450-1000 bp range. The shape and prominence of the trace is dependent on cell type, sublibrary size, and amount of DNA loaded into the TapeStation. Experimental conditions can also contribute to variability in the trace. Sublibrary with minor deviations can still produce high quality data.

Appendices

Appendix A: Sublibrary Generation Table

The Sublibrary Generation Table provides guidance for sample dilution in Section 1.5 to achieve a target sublibrary concentration.

Green text (top): Volume of cell suspension in μL (from step 1.5.10) to add to each sublibrary.

Purple text (bottom): Volume of Pre-Lysis Dilution Buffer in μL to add to each sublibrary.

Blue Shading: Serial dilution of cell stock is required to improve sublibrary cell count.

Red Shading: Insufficient cell stock concentration for target sublibrary cell count.

| Stock con (cells/ μL) | Target Sublibrary Cell Count (cells/sublibrary) | | | | | | | | | | | | | | |
|--------------------------------------|---|------|-------|-------|-------|--------|-------|--------|-------|--------|-------|--------|--------|--------|--------|
| | 200 | 500 | 1,000 | 2,000 | 3,000 | 4,000 | 5,000 | 6,000 | 7,000 | 8,000 | 9,000 | 10,000 | 11,000 | 12,000 | 12,500 |
| 50 | 4 | 10 | 20 | | | | | | | | | | | | |
| | 25 | 25 | 23 | | | | | | | | | | | | |
| 100 | 2 | 5 | 10 | 20 | | | | | | | | | | | |
| | 23 | 20 | 15 | 5 | | | | | | | | | | | |
| 200 | | 2.5 | 5 | 10 | 15 | 20 | 25 | | | | | | | | |
| | | 22.5 | 20 | 15 | 10 | 5 | 0 | | | | | | | | |
| 300 | | | 3.33 | 6.67 | 10 | 13.33 | 16.67 | 20 | 23.33 | | | | | | |
| | | | 21.67 | 18.33 | 15 | 11.67 | 8.33 | 5 | 1.67 | | | | | | |
| 400 | | | 2.5 | 5 | 7.5 | 10 | 12.5 | 15 | 17.5 | 20 | 22.5 | 25 | | | |
| | | | 22.5 | 20 | 17.5 | 15 | 12.5 | 10 | 7.5 | 5 | 2.5 | 0 | | | |
| 500 | | | 2 | 4 | 6 | 8 | 10 | 12 | 14 | 16 | 18 | 20 | 22 | 24 | 25 |
| | | | 23 | 21 | 19 | 17 | 15 | 13 | 11 | 9 | 7 | 5 | 3 | 1 | 0 |
| 600 | | | | 3.33 | 5 | 6.67 | 8.33 | 10 | 11.67 | 13.33 | 15 | 16.67 | 18.33 | 20 | 20.83 |
| | | | | 21.67 | 20 | 18.33 | 16.67 | 15 | 13.33 | 11.67 | 10 | 8.33 | 6.67 | 5 | 4.17 |
| 700 | | | | 2.86 | 4.29 | 5.71 | 7.14 | 8.57 | 10 | 11.43 | 12.86 | 14.29 | 15.71 | 17.14 | 17.86 |
| | | | | 22.14 | 20.71 | 19.29 | 17.86 | 16.43 | 15 | 13.57 | 12.14 | 10.71 | 9.29 | 7.86 | 7.14 |
| 800 | | | | 2.5 | 3.75 | 5 | 6.25 | 7.5 | 8.75 | 10 | 11.25 | 12.5 | 13.75 | 15 | 15.625 |
| | | | | 22.5 | 21.25 | 20 | 18.75 | 17.5 | 16.25 | 15 | 13.75 | 12.5 | 11.25 | 10 | 9.375 |
| 900 | | | | 2.22 | 3.33 | 4.44 | 5.56 | 6.67 | 7.78 | 8.89 | 10 | 11.11 | 12.22 | 13.33 | 13.89 |
| | | | | 22.78 | 21.67 | 20.56 | 19.44 | 18.33 | 17.22 | 16.11 | 15 | 13.89 | 12.78 | 11.67 | 11.11 |
| 1,000 | | | | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 12.5 |
| | | | | 23 | 22 | 21 | 20 | 19 | 18 | 17 | 16 | 15 | 14 | 13 | 12.5 |
| 1,100 | | | | | 2.73 | 3.64 | 4.55 | 5.45 | 6.36 | 7.27 | 8.18 | 9.09 | 10 | 10.91 | 11.36 |
| | | | | | 22.27 | 21.36 | 20.45 | 19.55 | 18.64 | 17.73 | 16.82 | 15.91 | 15 | 14.09 | 13.64 |
| 1,200 | | | | | 2.50 | 3.33 | 4.17 | 5 | 5.83 | 6.67 | 7.5 | 8.33 | 9.17 | 10 | 10.42 |
| | | | | | 22.50 | 21.67 | 20.83 | 20 | 19.17 | 18.33 | 17.5 | 16.67 | 15.83 | 15 | 14.58 |
| 1,300 | | | | | 2.31 | 3.08 | 3.85 | 4.62 | 5.38 | 6.15 | 6.92 | 7.69 | 8.46 | 9.23 | 9.62 |
| | | | | | 22.69 | 21.92 | 21.15 | 20.38 | 19.62 | 18.85 | 18.08 | 17.31 | 16.54 | 15.77 | 15.38 |
| 1,400 | | | | | 2.14 | 2.86 | 3.57 | 4.29 | 5 | 5.71 | 6.43 | 7.14 | 7.86 | 8.57 | 8.93 |
| | | | | | 22.86 | 22.14 | 21.43 | 20.71 | 20 | 19.29 | 18.57 | 17.86 | 17.14 | 16.43 | 16.07 |
| 1,500 | | | | | 2.00 | 2.67 | 3.33 | 4.00 | 4.67 | 5.33 | 6.00 | 6.67 | 7.33 | 8.00 | 8.33 |
| | | | | | 23.00 | 22.33 | 21.67 | 21.00 | 20.33 | 19.67 | 19.00 | 18.33 | 17.67 | 17.00 | 16.67 |
| 1,600 | | | | | 2.5 | 3.125 | 3.75 | 4.375 | 5 | 5.625 | 6.25 | 6.88 | 7.50 | 8.125 | 8.375 |
| | | | | | 22.5 | 21.875 | 21.25 | 20.625 | 20 | 19.375 | 18.75 | 18.13 | 17.50 | 17.19 | 16.875 |
| 1,700 | | | | | 2.35 | 2.94 | 3.53 | 4.12 | 4.71 | 5.29 | 5.88 | 6.47 | 7.06 | 7.65 | 7.95 |
| | | | | | 22.65 | 22.06 | 21.47 | 20.88 | 20.29 | 19.71 | 19.12 | 18.53 | 17.94 | 17.65 | 17.35 |
| 1,800 | | | | | 2.22 | 2.78 | 3.33 | 3.89 | 4.44 | 5 | 5.56 | 6.11 | 6.67 | 7.22 | 7.5 |
| | | | | | 22.78 | 22.22 | 21.67 | 21.11 | 20.56 | 20 | 19.44 | 18.89 | 18.33 | 18.06 | 17.75 |
| 1,900 | | | | | 2.11 | 2.63 | 3.16 | 3.68 | 4.21 | 4.74 | 5.26 | 5.79 | 6.32 | 6.84 | 7.14 |
| | | | | | 22.89 | 22.37 | 21.84 | 21.32 | 20.79 | 20.26 | 19.74 | 19.21 | 18.68 | 18.42 | 18.12 |
| 2,000 | | | | | 2 | 2.5 | 3.00 | 3.5 | 4.00 | 4.50 | 5 | 5.50 | 6.00 | 6.25 | 6.125 |
| | | | | | 23 | 22.5 | 22.00 | 21.5 | 21.00 | 20.50 | 20 | 19.50 | 19.00 | 18.75 | 18.5 |
| 2,100 | | | | | | 2.38 | 2.86 | 3.33 | 3.81 | 4.29 | 4.76 | 5.24 | 5.71 | 6.18 | 6.48 |
| | | | | | | 22.62 | 22.14 | 21.67 | 21.19 | 20.71 | 20.24 | 19.76 | 19.29 | 19.05 | 18.75 |
| 2,200 | | | | | | 2.27 | 2.73 | 3.18 | 3.64 | 4.09 | 4.55 | 5.00 | 5.45 | 5.88 | 6.18 |
| | | | | | | 22.73 | 22.27 | 21.82 | 21.36 | 20.91 | 20.45 | 20.00 | 19.55 | 19.32 | 19.02 |
| 2,300 | | | | | | 2.17 | 2.61 | 3.04 | 3.48 | 3.91 | 4.35 | 4.78 | 5.22 | 5.65 | 5.95 |
| | | | | | | 22.83 | 22.39 | 21.96 | 21.52 | 21.09 | 20.65 | 20.22 | 19.78 | 19.57 | 19.27 |

Appendix B: Sequencing Information

We recommend a minimum sequencing depth of 20,000 reads per cell for whole transcriptome libraries. However, ideal sequencing depth is dependent on the sample type and experimental goals. For example, if trying to capture rare cell types in a heterogenous population, deeper sequencing is required. Alternatively, analyzing a homogenous population, more shallow sequencing may be appropriate.

Sequencing libraries should be diluted and denatured according to the manufacturer’s instructions for the relevant sequencing instrument. If using the Illumina platform, we strongly recommend adding 5% PhiX for optimal sequencing quality.

Details of the final whole transcriptome sequencing library structure are below.




Figure 7: Whole transcriptome sequencing library structure.

Evercode sequencing libraries should be diluted and denatured according to the instruction for the relevant sequencing instrument. We strongly recommended adding 5% PhiX for optimal sequencing quality.

Details of the final TCR sequencing library structure are below.




Figure 8: TCR sequencing library structure.

Paired end sequencing is used for both Whole Transcriptome and TCR Library. Whole Transcriptome and TCR Libraries have different sequencing configurations as shown below.

| WHOLE TRANSCRIPTOME ILLUMINA RUN CONFIGURATION | | | |
|--|-----------------|--------------------|---|
| Run | Function | WT Library Cycles | Notes |
| Read 1 | cDNA Insert | 64 | The DNA product used for this sequencing run should have come from Section 3. |
| i7 Index (Index 1) | Barcode 4 / UDI | 8 | |
| i5 Index (Index 2) | Barcode 4 / UDI | 8 | |
| Read 2 | Barcodes 1-3 | 58 | |
| TCR ILLUMINA RUN CONFIGURATION | | | |
| Run | Function | TCR Library Cycles | Notes |
| Read 1 | CDR3 Insert | 242 | The DNA product used for this sequencing run should have come from Section 4. |
| i7 Index (Index 1) | Barcode 4 / UDI | 8 | |
| i5 Index (Index 2) | Barcode 4 / UDI | 8 | |
| Read 2 | Barcodes 1-3 | 58 | |

Appendix B1: Unique Dual Index (UDI) Plate - WT Sequences

Please refer to the following table to demultiplex the whole transcriptome libraries.

| Sublibrary Index ID | Well Position | i7 Forward Sequence | i5 Reverse Complement Sequence | i5 Forward Sequence |
|---------------------|---------------|---------------------|--------------------------------|---------------------|
| UDI_Plate_WT_1 | A1 | CAGATCAC | ATGTGAAG | CTTCACAT |
| UDI_Plate_WT_2 | B1 | ACTGATAG | GTCCAACC | GGTTGGAC |
| UDI_Plate_WT_3 | C1 | GATCAGTC | AGAGTCAA | TTGACTCT |
| UDI_Plate_WT_4 | D1 | CTTGTAAT | AGTTGGCT | AGCCAAC T |
| UDI_Plate_WT_5 | E1 | AGTCAAGA | ATAAGGCG | CGCCTTAT |
| UDI_Plate_WT_6 | F1 | CCGTCCTA | CCGTACAG | CTGTACGG |
| UDI_Plate_WT_7 | G1 | GTAGAGTA | CATTCATG | CATGAATG |
| UDI_Plate_WT_8 | H1 | GTCCGCCT | AGATACGG | CCGTATCT |
| UDI_Plate_WT_9 | A2 | GTGAAACT | TACAGACT | AGTCTGTA |
| UDI_Plate_WT_10 | B2 | TCATTCCT | AATGCCTG | CAGGCATT |
| UDI_Plate_WT_11 | C2 | GGTAGCAT | TGCTTGCC | GGCAAGCA |
| UDI_Plate_WT_12 | D2 | ACTTGATC | TTTGGGTG | CACCCAAA |
| UDI_Plate_WT_13 | E2 | ATGAGCAT | GAATCTGA | TCAGATTC |
| UDI_Plate_WT_14 | F2 | GCGCTATC | CGACTGGA | TCCAGTCG |
| UDI_Plate_WT_15 | G2 | TGACCAGT | ACATTGGC | GCCAATGT |
| UDI_Plate_WT_16 | H2 | TATAATCA | ACCACTGT | ACAGTGGT |
| UDI_Plate_WT_17 | A3 | CAAAAGTC | CGGTTGTT | AACAACCG |
| UDI_Plate_WT_18 | B3 | CGATGTCA | CATGAGGA | TCCTCATG |
| UDI_Plate_WT_19 | C3 | CTCAGAGT | TGGAGAGT | ACTCTCCA |
| UDI_Plate_WT_20 | D3 | TAATCGAC | TGACTTCG | CGAAGTCA |
| UDI_Plate_WT_21 | E3 | CATTTTCT | GGAAGGAT | ATCCTTCC |

| | | | | |
|-----------------|----|----------|-----------|----------|
| UDI_Plate_WT_22 | F3 | CTATACTC | TGTTTCGAG | CTCGAACA |
| UDI_Plate_WT_23 | G3 | CACTCACA | AAGGCTGA | TCAGCCTT |
| UDI_Plate_WT_24 | H3 | CTCGAACA | CTCGAGTG | CACTCGAG |
| UDI_Plate_WT_25 | A4 | CTCTATCG | ATCGGTGG | CCACCGAT |
| UDI_Plate_WT_26 | B4 | TCCTCATG | AGGTCTTG | CAAGACCT |
| UDI_Plate_WT_27 | C4 | AACAACCG | AGGAAGCG | CGCTTCCT |
| UDI_Plate_WT_28 | D4 | GCCAATGT | ACATGTGT | ACACATGT |
| UDI_Plate_WT_29 | E4 | TGGTTGTT | ATACAGTT | AACTGTAT |
| UDI_Plate_WT_30 | F4 | TCTGCTGT | ATCGCCTT | AAGGCGAT |
| UDI_Plate_WT_31 | G4 | TTGGAGGT | TTCGACGC | GCGTCGAA |
| UDI_Plate_WT_32 | H4 | TCGAGCGT | TGTCGTTC | GAACGACA |
| UDI_Plate_WT_33 | A5 | TGCGATCT | TCCATAGC | GCTATGGA |
| UDI_Plate_WT_34 | B5 | TTCCTGCT | TAAGTGTC | GACACTTA |
| UDI_Plate_WT_35 | C5 | TTCCATTG | CTGGCATA | TATGCCAG |
| UDI_Plate_WT_36 | D5 | TAACGCTG | CTGAGCCA | TGGCTCAG |
| UDI_Plate_WT_37 | E5 | TTGGTATG | CTCAATGA | TCATTGAG |
| UDI_Plate_WT_38 | F5 | TGAACTGG | CGCATACA | TGTATGCG |
| UDI_Plate_WT_39 | G5 | TCCAGTCG | CCGAAGTA | TACTTCGG |
| UDI_Plate_WT_40 | H5 | TGTATGCG | CCAGTTCA | TGAACTGG |
| UDI_Plate_WT_41 | A6 | TGGCTCAG | CAGCGTTA | TAACGCTG |
| UDI_Plate_WT_42 | B6 | TATGCCAG | CAATGGAA | TTCCATTG |
| UDI_Plate_WT_43 | C6 | GGTTGGAC | ATCCTGTA | TACAGGAT |
| UDI_Plate_WT_44 | D6 | GACACTTA | AGCAGGAA | TTCCTGCT |
| UDI_Plate_WT_45 | E6 | GAACGACA | ACGCTCGA | TCGAGCGT |
| UDI_Plate_WT_46 | F6 | AAGGCGAT | ACAGCAGA | TCTGCTGT |

| | | | | |
|-----------------|----|----------|----------|----------|
| UDI_Plate_WT_47 | G6 | ATGCTTGA | ACAAGCTA | TAGCTTGT |
| UDI_Plate_WT_48 | H6 | AGTATCTG | CATCAAGT | ACTTGATG |

Appendix B2: Unique Dual Index (UDI) Plate- EC Sequences

Please refer to the following table to demultiplex TCR sequencing libraries that have been sequenced together in the same run.

| Sublibrary Index ID | Well Position | i7 Forward Sequence | i5 Reverse Complement Sequence | i5 Forward Sequence |
|---------------------|---------------|---------------------|--------------------------------|---------------------|
| UDI_Plate_EC_1 | A1 | AACAGATC | TGCCATGA | TCATGGCA |
| UDI_Plate_EC_2 | B1 | CCTGTCTT | CTAACGAT | ATCGTTAG |
| UDI_Plate_EC_3 | C1 | ATATCGAG | TTCCGGAA | TTCCGGAA |
| UDI_Plate_EC_4 | D1 | TAGTCCGC | TCTCCATA | TATGGAGA |
| UDI_Plate_EC_5 | E1 | TGCTGTTA | AGTGACGT | ACGTCACT |
| UDI_Plate_EC_6 | F1 | TAGCGAAT | CTCTGGCA | TGCCAGAG |
| UDI_Plate_EC_7 | G1 | AGGACCGT | AGAAGATG | CATCTTCT |
| UDI_Plate_EC_8 | H1 | TTATCAGG | TCTAACAG | CTGTTAGA |
| UDI_Plate_EC_9 | A2 | CCTCGGAA | TACGTAGA | TCTACGTA |
| UDI_Plate_EC_10 | B2 | TTGATCGA | TAAGCTCT | AGAGCTTA |
| UDI_Plate_EC_11 | C2 | TCAACCTC | TTGGTCGA | TCGACCAA |
| UDI_Plate_EC_12 | D2 | CGGAATAA | CGGTAATA | TATTACCG |
| UDI_Plate_EC_13 | E2 | TATGAGAC | CACAGTGG | CCACTGTG |
| UDI_Plate_EC_14 | F2 | CCTACCAT | AGACGCGA | TCGCGTCT |
| UDI_Plate_EC_15 | G2 | CATAGGCC | AGCAATGG | CCATTGCT |
| UDI_Plate_EC_16 | H2 | CACGATAA | CCTCTTGA | TCAAGAGG |
| UDI_Plate_EC_17 | A3 | TTCTGGCG | TTGTCTAG | CTAGACAA |

| Sublibrary Index ID | Well Position | i7 Forward Sequence | i5 Reverse Complement Sequence | i5 Forward Sequence |
|---------------------|---------------|---------------------|--------------------------------|---------------------|
| UDI_Plate_EC_18 | B3 | TCGTA ACT | TGTTGCTA | TAGCAACA |
| UDI_Plate_EC_19 | C3 | TCCGTGGT | TACTGTCA | TGACAGTA |
| UDI_Plate_EC_20 | D3 | AAGACAAC | AAGCCATA | TATGGCTT |
| UDI_Plate_EC_21 | E3 | ATTGCTTC | CGGCACTT | AAGTGCCG |
| UDI_Plate_EC_22 | F3 | AGAGAAGT | TCATGGAT | ATCCATGA |
| UDI_Plate_EC_23 | G3 | TACGTCTT | TAATCAGG | CCTGATTA |
| UDI_Plate_EC_24 | H3 | TCCACGTT | CATCTGAG | CTCAGATG |
| UDI_Plate_EC_25 | A4 | CAATGAGT | CAGAATCG | CGATTCTG |
| UDI_Plate_EC_26 | B4 | CGTCAGTT | TTGATGCA | TGCATCAA |
| UDI_Plate_EC_27 | C4 | TTCACTCC | ACGTATGG | CCATACGT |
| UDI_Plate_EC_28 | D4 | CTACCTGA | TGCGGATA | TATCCGCA |
| UDI_Plate_EC_29 | E4 | AGGTGATT | TAACAGCA | TGCTGTTA |
| UDI_Plate_EC_30 | F4 | TCTCACAT | AATGTTCG | CGAACATT |
| UDI_Plate_EC_31 | G4 | TTAGTGAG | TTACGTGT | ACACGTAA |
| UDI_Plate_EC_32 | H4 | TGTTCACT | ATCTACTG | CAGTAGAT |
| UDI_Plate_EC_33 | A5 | CGAAGCCT | ATACCTCT | AGAGGTAT |
| UDI_Plate_EC_34 | B5 | TGGCCGTA | TACATCTG | CAGATGTA |
| UDI_Plate_EC_35 | C5 | CGCTTCAC | AAGGAGCA | TGCTCCTT |
| UDI_Plate_EC_36 | D5 | CGGAGAAC | CTCCTAGA | TCTAGGAG |
| UDI_Plate_EC_37 | E5 | AAGTTCAG | TCAGCCTA | TAGGCTGA |
| UDI_Plate_EC_38 | F5 | AAGCCTTC | TTCGCTCA | TGAGCGAA |
| UDI_Plate_EC_39 | G5 | ACACCTCA | CGGATTAA | TTAATCCG |
| UDI_Plate_EC_40 | H5 | TTCTCCTA | ATTAGAGG | CCTCTAAT |

| Sublibrary Index ID | Well Position | i7 Forward Sequence | i5 Reverse Complement Sequence | i5 Forward Sequence |
|---------------------|---------------|---------------------|--------------------------------|---------------------|
| UDI_Plate_EC_41 | A6 | TGACACGC | CGTGTGAA | TTCACACG |
| UDI_Plate_EC_42 | B6 | AAGCGCCT | TGTCACGG | CCGTGACA |
| UDI_Plate_EC_43 | C6 | TGATAACC | CGTAATCT | AGATTACG |
| UDI_Plate_EC_44 | D6 | AAGAGTGT | CCTTCTGG | CCAGAAGG |
| UDI_Plate_EC_45 | E6 | TGCTACCG | AATCGCTA | TAGCGATT |
| UDI_Plate_EC_46 | F6 | CCATCGTC | TGGACCAA | TTGGTCCA |
| UDI_Plate_EC_47 | G6 | CCACGAGA | TGTGCACT | AGTGCACA |
| UDI_Plate_EC_48 | H6 | ACGGTCAT | CTTATGGA | TCCATAAG |

Appendix C: Thermocycling Programs

Section 1: *In Situ* Cell Barcoding

| THAW ROUND 1 PLATE | | |
|--------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 10 min | 70°C | 26 µL |
| Step | Time | Temperature |
| 1 | 10 min | 25°C |
| 2 | Hold | 4°C |

| BARCODING ROUND 1 | | | |
|-------------------|-----------------|---------------|--------|
| Run Time | Lid Temperature | Sample Volume | |
| 40 min | 70°C | 40 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 10 min | 50°C | 1 |
| 2 | 12 s | 8°C | 3 |
| 3 | 45 s | 15°C | |
| 4 | 45 s | 20°C | |
| 5 | 30 s | 30°C | |
| 6 | 2 min | 42°C | |
| 7 | 3 min | 50°C | 1 |
| 8 | 5 min | 50°C | |
| 9 | Hold | 4°C | 1 |

| THAW ROUND 2 PLATE | | |
|--------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 10 min | 70°C | 10 µL |
| Step | Time | Temperature |
| 1 | 10 min | 25°C |
| 2 | Hold | 4°C |

| BARCODING ROUND 2 | | |
|-------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 50°C | 50 µL |
| Step | Time | Temperature |
| 1 | 15 min | 16°C |
| 2 | Hold | 4°C |

| ROUND 2 STOP | | |
|--------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 5 min | 50°C | 60 µL |
| Step | Time | Temperature |
| 1 | 5 min | 16°C |
| 2 | Hold | 4°C |

| THAW ROUND 3 PLATE | | |
|--------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 10 min | 70°C | 10 µL |
| Step | Time | Temperature |
| 1 | 10 min | 25°C |
| 2 | Hold | 4°C |

| BARCODING ROUND 3 | | |
|-------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 50°C | 60 µL |
| Step | Time | Temperature |
| 1 | 15 min | 16°C |
| 2 | Hold | 4°C |

| CELL LYSIS | | |
|------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 80°C | 55 µL |
| Step | Time | Temperature |
| 1 | 15 min | 65°C |
| 2 | Hold | 4°C |

Section 2: cDNA Capture and Amplification

| TEMPLATE SWITCH | | |
|-----------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 60 min | 70°C | 100 µL |
| Step | Time | Temperature |
| 1 | 60 min | 42°C |
| 2 | Hold | 4°C |

| cDNA AMPLIFICATION | | | |
|--------------------|-----------------|---------------|----------|
| Run Time | Lid Temperature | Sample Volume | |
| 50-70 min | 105°C | 100 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | 5 |
| 3 | 45 s | 65°C | |
| 4 | 3 min | 72°C | |
| 5 | 20 s | 98°C | |
| 6 | 20 s | 67°C | Variable |
| 7 | 3 min | 72°C | |
| 8 | 5 min | 72°C | 1 |
| 9 | Hold | 4°C | 1 |

Section 3: Sequencing Library Preparation

| FRAGMENTATION AND END PREP | | |
|----------------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 40 min | 70°C | 50 µL |
| Step | Time | Temperature |
| 1 | Hold | 4°C |
| 2 | 10 min | 32°C |
| 3 | 30 min | 65°C |
| 4 | Hold | 4°C |

| ADAPTER LIGATION | | |
|------------------|-----------------|---------------|
| Run Time | Lid Temperature | Sample Volume |
| 15 min | 30°C | 100 µL |
| Step | Time | Temperature |
| 1 | 15 min | 20°C |
| 2 | Hold | 4°C |

| INDEXING PCR | | | |
|--------------|-----------------|---------------|--|
| Run Time | Lid Temperature | Sample Volume | |
| ~30 min | 105°C | 50 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | Variable, see table in Section 3.5 |
| 3 | 20 s | 67°C | |
| 4 | 1 min | 72°C | |
| 5 | 5 min | 72°C | 1 |
| 6 | Hold | 4°C | 1 |

Section 4: TCR Sequencing Library Preparation

| TCR AMPLIFICATION 1 | | | |
|---------------------|-----------------|---------------|--------|
| Run Time | Lid Temperature | Sample Volume | |
| ~50 minutes | 105°C | 100 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | 20 |
| 3 | 20 s | 65°C | |
| 4 | 20 s | 72°C | |
| 5 | 5 min | 72°C | 1 |
| 6 | Hold | 4°C | 1 |

| TCR AMPLIFICATION 2 | | | |
|---------------------|-----------------|---------------|--------|
| Run Time | Lid Temperature | Sample Volume | |
| ~30 minutes | 105°C | 100 µL | |
| Step | Time | Temperature | Cycles |
| 1 | 3 min | 95°C | 1 |
| 2 | 20 s | 98°C | 9 |
| 3 | 20 s | 67°C | |
| 4 | 1 min | 72°C | |
| 5 | 5 min | 72°C | 1 |
| 6 | Hold | 4°C | 1 |

Appendix D: TCR-Specific Troubleshooting

Adjusting TCR Amplification Cycles

Section 4.1: TCR Amplification 1 and Section 4.3: TCR Amplification 2 requires 10 ng of input for optimal TCR amplification. However, PCR cycles can be adjusted to better suit different input quantities. See the tables below for recommended PCR cycles based on input.

| TCR AMPLIFICATION 1 (STEP 4.1.11) | |
|-----------------------------------|------------------------|
| Input (ng) | Total Number of Cycles |
| 10 | 20 |
| 5-9 | 21 |

| TCR AMPLIFICATION 2 (STEP 4.3.12) | |
|-----------------------------------|------------------------|
| Input (ng) | Total Number of Cycles |
| 10 | 9 |
| 5-9 | 10 |



Note: Having less than 5 nanograms of input material for TCR Amplification 1 or 2 might indicate poor quality data.

Appendix E: Revision History

| Version | Description | Date |
|---------|--|---------------|
| 1.0 | Updated for Evercode WT v3 chemistry | May 2024 |
| 1.1 | Appendix B: Updated sequencing information | June 2024 |
| 1.2 | Updated box configurations | August 2024 |
| 1.3 | Section 1.2.17: added a step to accommodate low Input fixation samples | November 2024 |



parsebiosciences.com

support@parsebiosciences.com

